
DK50 2V/M

- EN USER MANUAL
- DE BENUTZERHANDBUCH
- FR MANUEL DE L'UTILISATEUR
- RU РУКОВОДСТВО ПОЛЬЗОВАТЕЛЯ
- PL INSTRUKCJA OBSŁUGI
- SK NÁVOD NA POUŽITIE
- CS NÁVOD K POUŽITÍ

COMPRESSOR
KOMPRESSOR
COMPRESSEUR
КОМПРЕССОР
SPREŽARKA
KOMPRESOR

DK50 2V

EKOM spol. s r. o.
Priemyselná 5031/18
SK-921 01 Piešťany
Slovak Republic
tel.: +421 33 7967255
fax: +421 33 7967223

www.ekom.sk
email: ekom@ekom.sk

DATE OF LAST REVISION
DATUM DER LETZTEN ÜBERARBEITUNG
DATE DE LA DERNIÈRE RÉVISION
ДАТА ПОСЛЕДНЕЙ РЕДАКЦИИ
DATA OSTATNIEJ AKTUALIZACJI
DÁTUM POSLEDNEJ REVÍZIE
DATUM POSLEDNÍ REVIZE

09/2024

NP-DK50 2V-7_09-2024-AD
112000576-000

EN	<u>CONTENTS..... 5</u>
DE	<u>INHALT 62</u>
FR	<u>TABLES DES MATIÈRES 120</u>
RU	<u>ОГЛАВЛЕНИЕ 178</u>
PL	<u>SPIS TREŚCI 236</u>
SK	<u>OBSAH 294</u>
CS	<u>OBSAH 351</u>

CONTENTS

IMPORTANT INFORMATION	6
1. CONFORMITY WITH THE REQUIREMENTS OF EUROPEAN UNION DIRECTIVES	6
2. INTENDED USE	6
3. CONTRAINDICATIONS AND SIDE-EFFECTS	6
4. SYMBOLS.....	6
5. WARNINGS	7
6. STORAGE AND TRANSPORT	9
PRODUCT DESCRIPTION	10
7. VARIANTS	10
8. ACCESSORIES	10
9. PRODUCT FUNCTION	11
TECHNICAL DATA	16
INSTALLATION	25
10. INSTALLATION CONDITIONS	25
11. COMPRESSOR ASSEMBLY	25
12. PNEUMATIC CONNECTION	27
13. ELECTRICAL CONNECTION	30
14. COMPRESSOR PLACEMENT IN THE CABINET	31
15. COMMISSIONING	33
16. PNEUMATIC AND ELECTRICAL DIAGRAMS	34
OPERATION	40
17. SWITCHING ON THE COMPRESSOR.....	41
18. SWITCHING OFF THE COPRESSOR	41
AD DRYER	42
19. PRINCIPLE OF OPERATION	42
20. USER INTERFACE/SETTINGS	43
21. MAIN SCREENS.....	45
22. INDICATION OF MAINTENANCE AND ALARM STATES.....	47
PRODUCT MAINTENANCE	50
23. PRODUCT MAINTENANCE.....	50
24. LONG-TERM SHUTDOWN.....	58
25. DISPOSAL OF DEVICE	58
TROUBLESHOOTING	59
26. REPAIR SERVICE	61
ANNEX	408
27. INSTALLATION RECORD	408

IMPORTANT INFORMATION

1. CONFORMITY WITH THE REQUIREMENTS OF EUROPEAN UNION DIRECTIVES

This product conforms to the requirements of the Regulation (EU) on medical devices (MDR 2017/745) and is safe for the intended use if all safety instructions are followed.

2. INTENDED USE

The compressor is used as a source of clean, oil-free compressed air to power active medical devices where the parameters and properties of the compressed air are suitable for the specific application.

Compressed air supplied by the compressor is unsuitable for use with artificial lung ventilation devices

Any other use of the product beyond its intended use is considered an incorrect use. The manufacturer is not liable for any damages or injuries resulting from the incorrect use.

3. CONTRAINDICATIONS AND SIDE-EFFECTS

There are no contraindications or side-effects known.

4. SYMBOLS

The following symbols and marks are used in the User manual, on the device and its packaging:

General warning

Warning – risk of electric shock

Warning - compressor is controlled automatically

Warning – hot surface

General caution

Read the operating instructions

Refer to instruction manual

CE – marking



Medical device

Serial number

Article number

	Unique Device Identifier
	Swiss Authorised Representative
	Swiss Importer
	Protecting earthing
	Fuse
	Package handling label – fragile
	Package handling label – this side up
	Package handling label – keep dry
	Package handling label – temperature limits
	Package handling label – limited stacking
	Package label – recyclable material
	Manufacturer

5. WARNINGS

The product is designed and manufactured to be safe for the user and the surrounding environment when used in the defined manner. Keep the following warnings in mind.

5.1. General warnings

READ THE USER MANUAL CAREFULLY BEFORE USING THE DEVICE AND KEEP IT FOR FUTURE USE!

- The user manual supports the correct installation, operation and maintenance of the product. Careful review of this manual will provide the information necessary for the correct operation of the product for its intended use.
- Original packaging should be retained for the return of the device. Only the original packaging ensures protection of the device during transport. If it is necessary to return the product during the warranty period, the manufacturer is not liable for damages caused by improper packaging.
- This warranty does not cover damages originating from the use of accessories other than those specified or recommended by the manufacturer.
- The manufacturer only guarantees the safety, reliability and function of the device if:
 - installation, new settings, changes, expansion, and repairs are performed by the

manufacturer or an organisation authorised by the manufacturer

- the product is used pursuant to the user manual
- The user manual corresponds to the configuration of the product and its compliance with the applicable safety and technical standards at the time of printing. The manufacturer reserves all rights for the protection of its configuration, methods and names.
- Translation of the user manual is performed in accordance with the best available knowledge. The Slovak version is to be used in the event of any uncertainties.
- This user manual is the original instructions. Translation is performed in accordance with the best available knowledge.

5.2. General safety warnings

The manufacturer designed and manufactured the product to minimise all risks when used correctly for the intended use. The manufacturer considers it its obligation to lay down the following general safety precautions.

- Use and operation of the product must comply with all laws and local regulations valid in the place of use. The operator and user are responsible for following all the appropriate regulations in the interests of performing work safely.
- Only the use of original parts guarantees the safety of operating personnel and the flawless operation of the product itself. Only the accessories and parts mentioned in the technical documentation or expressly approved by the manufacturer should be used.
- The operator must ensure that the device is functioning correctly and safely before every use.
- The user must be familiar with the operation of the device.
- Do not use the product in environments with a risk of explosion.
- The user must inform the supplier immediately if any problem directly related to the operation of the device occurs.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

5.3. Safety warnings on protection from electric current

- The device must only be connected to a properly installed, earthed socket.
- Before the product is plugged in, ensure that the mains voltage and frequency stated on the product are in compliance with the values of the mains.
- Prior to putting the device into operation check for any damage to the connected pneumatic lines and electrical wiring. Replace damaged pneumatic lines and electrical wirings immediately.
- Immediately disconnect the product from the mains (remove power cord from the socket) in hazardous situations or in the case of a technical malfunction.
- During all repairs and maintenance, ensure that:
 - the mains plug is removed from the power socket
 - pressure is vented from the air tank and pipes
- The product shall only be installed by a qualified technician.

6. STORAGE AND TRANSPORT

The manufacturer ships the compressor in a transport packaging. This protects the device from damage during transport.

Risk of damage to pneumatic parts.

The compressor must be transported only when all air has been vented. Before moving or transporting the compressor, release all the air pressure from the tank and pressure hoses and drain condensate from the tank and from the condensate separator on the dryer.

Original packaging must be kept for eventual return of the device. If possible, always use the original compressor packaging for optimal protection of the product. If it is necessary to return the product within the guarantee period, the manufacturer does not guarantee for damages caused by incorrect packaging of the product.

Transport the compressor in an upright position, always secured by transport fixing elements.

Protect the compressor from high humidity, dirt and extreme temperatures during transport and storage. Do not store together with volatile chemical substances.

If it is not possible to keep the original packaging, dispose of the packaging in an environmentally friendly way. Transport cardboard can be recycled with old paper.

It is not permitted to store and transport the device outside the defined conditions, see below.

6.1. Ambient conditions

Products can be stored in rooms and means of transport that are free from any traces of volatile chemical substances under the following climatic conditions:


Temperature	from -25°C o $+55^{\circ}\text{C}$, in 24h to $+70^{\circ}\text{C}$
Relative humidity	max. 90% (without condensation/ non-condensing)

PRODUCT DESCRIPTION


7. VARIANTS

The compressor is manufactured according to its intended application in the following variants:

DK50 2V/M	compressor on base with air dryer
DK50 2V S/M	compressor in cabinet with air dryer

DK50 2V/M

DK50 2V S/M

8. ACCESSORIES

Accessories are not included in a standard order, they must be ordered separately.

Set of compressed air outlet filters

The compressor may be equipped with a set of filters if specified. The filter set may be equipped with an air pressure regulator.

Where a different level of air filtration is required, this requirement must be agreed upon with the supplier and specified in the order.

Type	Use	Level of filtration (μm)	Article number
FS 30F		1	
FS 30M	DK50 2V/M	1+0.1	to order
FS 30S		1+0.01	
FS 30A		1+0.01+A**	

**) A – submicrofilter with activated carbon

Filter set regulator assembly

The compressor may be equipped with a pressure regulator of the compressed air outlet if specified. The regulator must be selected according to the application to the filter set, or separately. The regulator shall ensure constant pressure at the outlet.

Type	Use	Article number
Regulator complete (to the filter assembly)	DK50 2V/M	to order
REG11 (separate)		447000001-077

Filter set brackets

A suitable bracket must be ordered for every filter set.

Type	Use	Article number
Compressor-mounted bracket	DK50 2V/M	to order

Equipotential bonding socket

The socket provides protective bonding.

Type	Use	Name	Article no.
POAG-KBT6-EC	DK50 2VS/M	Mains plug connector	033400075-000
FLEXI-S/POAG-HK6		Conductor (1 m)	034110083-000

9. PRODUCT FUNCTION

9.1. Compressor with adsorption dryer

(Fig. 1) The air pump (1) draws in air through the inlet filter (8) and supplies compressed air through the connecting hose into the adsorption dryer (38). From the inlet to the dryer module, the air is first cooled in the integrated cooler (26) and then moves through the condensate separator (32), entering the active chamber with adsorbent (27), where the air is dried. A portion of the air is fed into the second, regenerated chamber, where this air is used to remove moisture from the adsorbent, after which it is vented to the atmosphere. The active chamber switches on a cyclical basis. The dry air then passes through the non-return valve (3) into the air tank (2). Air drying efficiency (the pressure dew point) is determined by the dew point sensor located in the air tank.

It is not necessary to drain the air tank.

9.2. Compressor in cabinet

(Fig. 3) The soundproof cabinet provides compact covering of the compressor, dampens noise efficiently, providing sufficient exchange of cooling air. Its design makes it suitable for use in rooms where low noise levels are required. The fan under the compressor pump provides cooling of the compressor, it is in operation at the same time with compressor motor or after switching on the temperature switch at temperature above 40 °C. After cooling the cabinet area approximately under 32°C, the fans turn off automatically.

The dryer display (34) is located on the front panel of the cabinet.

Risk of compressor overheating.

Make sure that there are no obstacles at the cooling air inlet into the cabinet (around the bottom part of the cabinet) and at the hot air outlet on the top back side of the cabinet.

If the compressor is placed on a soft floor, e.g. carpet, create space between the base and the floor or the cabinet and the floor, e.g. underlay the footings with hard pads to ensure sufficient cooling of the compressor.

Description to Figures 1 - 3

1	Air pump	22	Handle
2	Air tank	23	Casters
3	Non-return valve	24	Pressure gauge hose
4	Pressure switch	25	Outlet module
5	Safety valve	26	Cooling module
6	Pressure gauge	27	Dryer chamber
7	Drain valve	28	Inlet valve module
8	Inlet filter	29	Compressed air inlet
9	Compressor fan	30	Regeneration solenoid valve
10	Condensate collection vessel	31	Dryer fan
11	Magnetic holder	32	Condensate separator
12	Electric panel	33	Condensate drain solenoid valve
13	Cabinet fan	34	*Display
14	Door pin	35	Air outlet
15	Rectifying screw	36	**Bushing
16	Lock	37	Outlet from condensate separator and regeneration valves
17	Supporting bracket	38	Adsorption dryer
18	Wall stop	39	Cabinet connector with cord
19	Switch	40	Circuit breaker switch
20	Cabinet pressure gauge	41	*** Dew point sensor
21	Door hinge		

* Position 34 applies to the DK50 2V/M model without the cabinet; position 34 in the case of the DK50 2VS/M with the cabinet is a simple cover.

** For the cabinet model, the cable connected to the display (36) in the front panel of the cabinet is routed through the universal bushing.

*** Applies to model with dew point sensor.

Fig. 1: DK50 2V/M – Compressor with adsorption dryer




Fig. 2: Adsorption dryer**ADM140/25**

- * Position 34 applies to the DK50 2V/M model without the cabinet; position 34 in the case of the DK50 2VS/M with the cabinet is a simple cover.
- ** For the cabinet model, the cable connected to the display (36) in the front panel of the cabinet is routed through the universal bushing.

Fig. 3: Cabinet

TECHNICAL DATA

Compressors are designed to operate in dry, ventilated and indoor dust-free rooms with the following climatic conditions:

Temperature from +5°C to +30°C

Relative humidity max. 50%


Working pressure 6 – 8 bar		DK50 2V/M	DK50 2V S/M
Nominal voltage Frequency ^{a)}	V, Hz	230, 50	230, 50
Capacity at 6 bar (FAD)	l/min	120	120
Working pressure ^{b)}	bar	6.0 – 8.0	6.0 – 8.0
Rated current	A	8.8	9.0
Motor power	kW	1.2	1.2
Air tank volume	l	25	25
Air quality – filtration	µm	-	-
Maximum operating pressure of safety valve	bar	11.0	11.0
Noise level at 5 bar (L _{pA})	dB	≤69.0	≤55.0
Operating mode	%	S1-100	S1-100
PDP drying at 7 bar	°C	≤ -20	≤ -20
Dryer low pressure alarm	bar	< 3	< 3
Dryer high pressure alarm	bar	> 9	> 9
Dryer high dew point alarm	°C	> -20	> -20
Time to fill air tank from 0 to 7 bar	s	90	90
Dimensions (net) w x d x h	mm	460x620x720	560x720x880
Net weight ^{c)}	kg	69	110
Classification under EN 60601-1	Class I.		

^{a)} State the compressor variant in the order

^{b)} For other range of pressure consult with the supplier

^{c)} Weight value is only informative data and applies only to a product without any accessories

Dependence of compressor output on working pressure


Working pressure 6 – 8 bar		DK50 2V/M	DK50 2V S/M
Nominal voltage, Frequency ^{a)}	V, Hz	230, 50	230, 50
Capacity at 6 bar (FAD)	l/min	115	115
Working pressure ^{b)}	bar	6.0 – 8.0	6.0 – 8.0
Rated current	A	8.8	9.0
Motor power	kW	1.2	1.2
Air tank volume	l	25	25
Air quality – filtration	µm	-	-
Maximum operating pressure of safety valve	bar	11.0	11.0
Noise level at 5 bar (L _{pA})	dB	≤69.0	≤55.0
Operating mode	%	S1-100	S1-100
PDP drying at 7 bar	°C	≤ -40	≤ -40
Dryer low pressure alarm	bar	< 3	< 3
Dryer high pressure alarm	bar	> 9	> 9
Dryer high dew point alarm	°C	> -40	> -40
Time to fill air tank from 0 to 7 bar	s	100	100
Dimensions (net) w x d x h	mm	460x620x720	560x720x880
Net weight ^{c)}	kg	69	110
Classification under EN 60601-1		Class I.	

^{a)} State the compressor variant in the order

^{b)} For other range of pressure consult with the supplier

^{c)} Weight value is only informative data and applies only to a product without any accessories

Dependence of compressor output on working pressure


Working pressure 6 – 8 bar		DK50 2V/M	DK50 2V S/M
Nominal voltage, Frequency ^{a)}	V, Hz	230, 50	230, 50
Capacity at 6 bar (FAD)	l/min	100	100
Working pressure ^{b)}	bar	6.0 – 8.0	6.0 – 8.0
Rated current	A	8.8	9.0
Motor power	kW	1.2	1.2
Air tank volume	l	25	25
Air quality – filtration	µm	-	-
Maximum operating pressure of safety valve	bar	11.0	11.0
Noise level at 5 bar (L _{pA})	dB	≤69.0	≤55.0
Operating mode	%	S1-100	S1-100
PDP drying at 7 bar	°C	≤ -70	≤ -70
Dryer low pressure alarm	bar	< 3	< 3
Dryer high pressure alarm	bar	> 9	> 9
Dryer high dew point alarm	°C	-68	-68
Time to fill air tank from 0 to 7 bar	s	110	110
Dimensions (net) w x d x h	mm	460x620x720	560x720x880
Net weight ^{c)}	kg	69	110
Classification under EN 60601-1	Class I.		

^{a)} State the compressor variant in the order

^{b)} For other range of pressure consult with the supplier

^{c)} Weight value is only informative data and applies only to a product without any accessories

Dependence of compressor output on working pressure


Working pressure 8 – 10 bar		DK50 2V/M	DK50 2V S/M
Nominal voltage, Frequency ^{a)}	V, Hz	230, 50	230, 50
Capacity at 8 bar (FAD)	l/min	100	100
Working pressure ^{b)}	bar	8.0 – 10.0	8.0 – 10.0
Rated current	A	8.6	8.8
Motor power	kW	1.2	1.2
Air tank volume	l	25	25
Air quality – filtration	µm	-	-
Maximum operating pressure of safety valve	bar	11.0	11.0
Noise level at 5 bar (L _{pA})	dB	≤69.0	≤55.0
Operating mode	%	S1-100	S1-100
PDP drying at 7 bar	°C	≤ -20	≤ -20
Dryer low pressure alarm	bar	< 3	< 3
Dryer high pressure alarm	bar	> 11	> 11
Dryer high dew point alarm	°C	> -20	> -20
Time to fill air tank from 0 to 7 bar	s	115	115
Dimensions (net) w x d x h	mm	460x620x720	560x720x880
Net weight ^{c)}	kg	71	112
Classification under EN 60601-1		Class I.	

^{a)} State the compressor variant in the order

^{b)} For other range of pressure consult with the supplier

^{c)} Weight value is only informative data and applies only to a product without any accessories

Dependence of compressor output on working pressure


Working pressure 8 – 10 bar		DK50 2V/M	DK50 2V S/M
Nominal voltage, Frequency ^{a)}	V, Hz	230, 50	230, 50
Capacity at 8 bar (FAD)	l/min	95	95
Working pressure ^{b)}	bar	8.0 – 10.0	8.0 – 10.0
Rated current	A	8.6	8.8
Motor power	kW	1.2	1.2
Air tank volume	l	25	25
Air quality – filtration	µm	-	-
Maximum operating pressure of safety valve	bar	11.0	11.0
Noise level at 5 bar (L _{pA})	dB	≤69.0	≤55.0
Operating mode	%	S1-100	S1-100
PDP drying at 7 bar	°C	≤ -40	≤ -40
Dryer low pressure alarm	bar	< 3	< 3
Dryer high pressure alarm	bar	> 11	> 11
Dryer high dew point alarm	°C	> -40	> -40
Time to fill air tank from 0 to 7 bar	s	120	120
Dimensions (net) w x d x h	mm	460x620x720	560x720x880
Net weight ^{c)}	kg	71	112
Classification under EN 60601-1	Class I.		

^{a)} State the compressor variant in the order

^{b)} For other range of pressure consult with the supplier

^{c)} Weight value is only informative data and applies only to a product without any accessories

Dependence of compressor output on working pressure

FAD correction of capacity for altitude

Capacity given in the form of FAD („Free Air Delivery“) applies to the following conditions:

Altitude	0 m.n.m.	Temperature	20°C
Atmospheric pressure	101325 Pa	Relative humidity	0%

To calculate FAD compressor capacity in dependence on altitude, it is necessary to apply correction factor according to the following table:

Altitude [m.n.m.]	0 - 1500	1501 - 2500	2501 - 3500	3501 - 4500
FAD correction factor	1	0,8	0,71	0,60

Electromagnetic compatibility declaration

Medical equipment needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided below.

Guidance and manufacturer's declaration - electromagnetic emissions		
Pursuant to IEC 60601-1-2:2014/AMD1:2020 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances		
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The equipment uses RF energy only for its internal functions. Therefore, the RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The equipment is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	The equipment is not likely to cause any flicker, as the current flow is approx. constant after the start up.	

Guidance and manufacturer's declaration - electromagnetic immunity			
Pursuant to IEC 60601-1-2:2014/AMD1:2020 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances			
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV 100 kHz repetition frequency Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV L-N ±2 kV L-PE; N-PE Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruption, and voltage variations on power supply input lines IEC 60601-4-11	$U_T=0\%$, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°) $U_T=0\%$, 1 cycle $U_T=70\%$ 25/30 cycles (at 0°) $U_T=0\%$, 250/300 cycles	$U_T=>95\%$, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°) $U_T=>95\%$, 1 cycle $U_T=70\%$ (30% dipp U_T), 25(50Hz)/30(60Hz) cycles (at 0°) $U_T=>95\%$, 250(50Hz)/300(60Hz) cycle	Mains power quality should be that of a typical commercial or hospital environment. The device stops and restarts automatically at each dip. In this case does not occur unacceptable pressure drop.
Power frequency (50/60 Hz) IEC 61000-4-8	N/A	N/A	The equipment does not contain magnetically sensitive components or circuits.
Radiated fields in close proximity per IEC 61000-4-39	N/A	N/A	The equipment does not contain magnetically sensitive components or circuits.
NOTE: U_T is the A.C. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic immunity			
Pursuant to IEC 60601-1-2:2014/AMD1:2020 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances			
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d=1,2\sqrt{P}$ $d=1,2\sqrt{P}, 80 \text{ MHz to } 800 \text{ MHz}$ $d=2,3\sqrt{P}, 800 \text{ MHz to } 2,7 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range ^b. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz	3 V/m	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment.			
^b Over the frequency range 150 kHz to 80 kHz, field strengths should be less than 3 V/m.			

INSTALLATION**Risk of incorrect installation.**

Only a qualified professional can install the compressor and place it into operation for the first time. His obligation is to train the operating personnel on the use and maintenance of the device. He shall confirm installation and training of operators by an entry into the installation record. (See warranty card)

10. INSTALLATION CONDITIONS

- The compressor may only be installed and operated in dry, well-ventilated and dust-free rooms with conditions stated in Technical data

**Risk of damage to the device.**

Do not operate the device in outdoor environment or in humid or wet rooms.

**Explosion risk.**

Do not install the device in any rooms, in which explosive gases, dusts or flammable liquids may be present.

- The compressor must be installed in such a way, that it is accessible at all times for operation and maintenance and label on the device is accessible.
- The compressor must stand on a flat, sufficiently stable base (be aware of the compressor weight, see Technical data).

You may notice a “new product” odour when you first place the product into service (for a short period of time). This odour is temporary and does not impede the normal use of the product. Ensure the space is properly ventilated after installation.

11. COMPRESSOR ASSEMBLY**11.1. Manipulation and fixation removal**

- Unpack the compressor from the packaging.
- For compressors with a cabinet, open the front door and remove the connecting strip (17) in the front lower part. Remove the cabinet from the compressor.
- Put the compressor at the place of operation (Fig. 4).



Use handling straps and the handle on the compressor to move the compressor as needed. Do not use other parts of the compressor (the air pump, cooler, etc.) for gripping.

The number of persons required to handle the equipment must be matched to the weight of the device.

Fig. 4: Manipulation with the compressor using handling straps

Securing of handling straps after putting compressor at the place of operation


- Remove transport fixing elements from the pumps (Fig. 5).

Prior to installation, ensure that the compressor is free of all transport packaging and stabilizers to avoid any risk of damage to the product.

Remove fixing elements of the pumps once the compressor has been assembled and balanced at the final installation place.

Fig. 5: Releasing the pump

12. PNEUMATIC CONNECTION

12.1. Connecting to the compressed air outlet

- Connect the pressure hose (2) to the compressed air outlet (1) of the compressor (Fig. 6).
- Secure the pressure hose in the clamp (3) (Fig. 6).
- Route the pressure hose through the opening on the left side in the back wall of the cabinet (Fig. 8) on compressors in cabinet.


Fig. 6: Connecting to the compressed air outlet

12.2. Connecting the cabinet pressure gauge to the compressor

- Remove the plug (1) from the threads (2) on the pneumatic block of the compressor.
- Connect the cabinet pressure gauge hose to the threaded fitting.


Fig. 7: Connecting the cabinet pressure gauge to the compressor

12.3. Connecting the condensate vessel


- Connect the condensate drain hose to the condensate bottle on compressors in cabinet. (Fig. 8)
- Put the plug in free inlet of the condensate collection vessel.
- Route the hose through the opening in the rear wall of the cabinet for compressors with dryers. Connect the condensate collection tank using the hose to the condensate outlet on the dryer.


Injury hazard.

Do not connect the condensate drain directly to a waste drain! Passers-by may be injured!

Fig. 8: Connection of condensate collection vessel


Risk of damage to pneumatic components.

Air hoses must not be broken.

12.4. Floor installation

Connect the compressor using the hoses prepared in advance in the floor, as per the installation plan.

Fig. 9: Floor installation

Description to Fig. 9:

- | | |
|---------------------------|-----------------------------------|
| 1. Outline of the cabinet | 4. Compressed air connection G3/8 |
| 2. Outline of the base | 5. Power supply input |
| 3. Front section - door | 230V/50(60) Hz |
| | 3G1.0x4000 |

- All dimensions are in millimetres
- Distance from wall min. 100 mm

13. ELECTRICAL CONNECTION

- The product is delivered with a cord equipped with a plug with earthing pin.
- Keep the socket easily accessible to ensure that the device can be safely disconnected from the mains.
- The corresponding circuit must be protected in the power distribution max. 16 A.

Risk of electric shock.

it is necessary to follow all local electrotechnical regulations. The mains voltage and frequency must comply with the data stated on the device label.

Risk of fire and electric shock.

Electrical cable must not be in contact with hot compressor components.

Risk of fire and electric shock.

Electrical cord must not be broken.

DK50 2V S/M compressor cabinet, variant 6-8 bar, is equipped with a cooling fan and a switch. The cabinet must be connected using the cord with a connector to a connector at the compressor electrical box.

- Connect the pin for equipotential bonding \varnothing 6 mm (1) to the electrical circuit using the defined method under valid electrical engineering regulations.
- The socket for equipotential bonding (2) is an accessory and is not provided in the basic package.

Fig. 10: Equipotential bonding


13.1. Connecting a compressor not installed in a cabinet

- Insert the mains plug into a rated mains socket.
- The compressor is ready for operation.

13.2. Connecting a compressor installed in a cabinet

- Route the mains plug through the opening in the rear wall of the cabinet for cabinet-mounted compressors.
- Connect the cabinet electrically to the compressor by inserting the provided power cord with connector (1) into an outlet. (Fig. 11)
- Disconnect the cabinet electrically from the compressor by pulling the connector (1) from the outlet while the latch is released.
- Condition the display at the front door of the cabinet and then secure the cord to the display in the clamps. (Fig. 16)

Fig. 11: Connecting a compressor installed in a cabinet

14. COMPRESSOR PLACEMENT IN THE CABINET

- Position the compressor within a distance of at least 500 mm from the cabinet to facilitate the movement of the hoses and the power cord within the cabinet (Fig. 12).


Fig. 12: Positioning of the compressor for connections

- Make the pneumatic connections per Chapter 12.1.
- Connect the condensate collection tank per Chapter 12.3.
- Route the pressure hose, the condensate drain hose, and the power cord through the opening in the rear wall of the cabinet.
- Position the compressor at such a distance from the cabinet so the WINSTA connector can make the electrical connection between the compressor and the cabinet (Fig. 11, Fig. 13).


Fig. 13: Positioning the compressor for electrical connections

- Electrically connect the compressor per Chapter 13.
- Slide the compressor into the cabinet and secure with the connecting strip (17).

A compressor with air dryer must not be inserted to the cabinet to the stop as there is a risk of permanent damage to the dryer.

- Check to ensure the compressor is fully inserted into the cabinet using the control dimension (Fig. 14).


Fig. 14: Positioning of the DK50 2V/M compressor in the cabinet

- Connect the cabinet pressure gauge to the compressor per Chapter 12.2.
- Connect the display in the front door of the cabinet (Fig. 15), and secure the cord to the display in the clamps (Fig. 16).

Fig. 15: Connecting the display connector


Fig. 16: Position of the clamps in the cabinet

15. COMMISSIONING

- Check if all fixing elements used during transport have been removed.
- Check correct connection of compressed air (see chapter 12).
- Check correct connection to the mains (see chapter 13)
- Check position of the circuit breaker switch, it must be in position „I“. If it is in position „0“, turn the switch (2) to position „I“ (Fig. 17)
- On compressor in the cabinet, switch the switch (5) on the front side of the cabinet to position „I“, the green light indicates the device status in operation. (Fig. 17).

- Check connection of the cabinet manometer hose to the pneumatic block of the compressor.
- Check the connection of the condensate collection vessel. (see chap. 12.3)
- Check the connection of the display at the front door of the cabinet.
- Set the language, time and drying mode (see chap. 20.2) on the compressor display.


The compressor is not equipped with a backup power supply.

16. PNEUMATIC AND ELECTRICAL DIAGRAMS

16.1. Pneumatic diagram

DK50 2V/M, DK50 2V S/M

Description to pneumatic diagrams:


- | | |
|------------------------|--------------------------------------|
| 1 Inlet filter | 10 Cooler |
| 2 Air pump | 11 Cooler fan |
| 3 Fan | 12 Air tank |
| 4 Safety valve | 13 Condensate drain valve |
| 5 Non-return valve | 14 Outlet valve |
| 6 Condensate separator | 15 Solenoid valve – regeneration |
| 7 Dryer | 16 Solenoid valve – condensate drain |
| 8 Pressure switch | 17 Condensate collection vessel |
| 9 Pressure gauge | |

16.2. Electrical diagram

DK50 2V/M, DK50 2V S/M

1/N//PE 230 V, 50 Hz

ELEKTRIC OBJECT OF 1st. CAT.

Note:


- Only connect the jumper for compressors without cabinet (see chap. 23.10)

AD140

1/N/PE ~ 230V, 50Hz / 230V, 60Hz

~ 120V, 60Hz

ELEKTRIC OBJECT OF 1st. CAT.


**Note:**

- The dew point sensor (A1_1) is not included with the dryer; please contact your supplier for more information.

AD280

1/N/PE ~ 230V, 50Hz / 230V, 60Hz
 ~ 120V, 60Hz

ELEKTRIC OBJECT OF 1st. CAT.

Note:


- The dew point sensor (A1_1) is not included with the dryer; please contact your supplier for more information.

AD500

1/N/PE ~ 230V, 50Hz / 230V, 60Hz

~ 120V, 60Hz

ELEKTRIC OBJECT OF 1st. CAT.


**Note:**

- The dew point sensor (A1_1) is not included with the dryer; please contact your supplier for more information.

Compressor cabinet

1/N/E 230 V, 50 Hz

ELEKTRIC OBJECT OF 1st. CAT.

Description to electrical diagrams:

M1	Compressor motor	F1	Circuit breaker switch
E1	Compressor fan	S10	Switch
E7-E9	Dryer fan	P11	Indicator
E10	Cabinet fan	RSA 10	Controller
B2	Temperature switch	M11, M12	Solenoid valve - regeneration
X10,X11	Connector	M4	Relief solenoid valve
B1	Pressure switch	X2, X3	Socket
X1	Terminal box	A1_1	Dew point sensor connector

OPERATION

THE EQUIPMENT MAY ONLY BE OPERATED BY THE TRAINED STAFF!

Risk of burn or fire.

During compressor operation, the pump parts may heat to temperatures dangerous for contact with persons or materials.

Warning – compressor is controlled automatically.

Automatic start-up. When pressure in the pressure tank drops to switch-on pressure, the compressor automatically switches on. The compressor automatically switches off, when pressure in the air tank reaches the switch-off pressure.

Hazard of damage to pneumatic components.

Working pressures of the pressure switch set by the manufacturer must not be changed. If the compressor operates at higher working pressure than the switch-on pressure, it means that the air consumption is too high (see Chapter Malfunctions)

Hazard of damage to dryer.

When the dryer is operated at ambient temperature that is higher than maximum operating temperature stated in Technical data, the dryer can be damaged.

Required drying performance can only be achieved when the defined operating conditions are followed.

Drying performance decreases and the achieved dew point drops if the dryer is operated at any pressure below the minimum working pressure.

During longer operation of the compressor, the temperature in the surrounding of the compressor increases above 40 °C and the cooling fan switches on automatically. After cooling the space below 32 °C, the fan switches off.

17. SWITCHING ON THE COMPRESSOR


Start the compressor (without a cabinet) at the pressure switch (1) by turning the switch (2) to position "I." This starts the compressor and fills the tank to the switching off pressure, which then shuts off the compressor.

Start the compressor (with a cabinet) at the switch (5) on the front side of the cabinet. This starts the compressor and fills the tank to the switching off pressure, which then shuts off the compressor.

The compressor operates in automatic mode, switched on and off by the pressure switch, depending on compressed air consumption.

Check the values of the switch-on and switch-off pressure using the pressure gauge (3). A tolerance of $\pm 10\%$ is acceptable. The pressure in the air tank must not exceed the permitted operating pressure.

Fig. 17: Switching on the compressor

- **Compressor** – when first started and placed into service, the compressor fills the air tank until the switching pressure and the pressure switch turn off the compressor. The compressor operates in automatic mode, switched on and off by the pressure switch, depending on compressed air usage.
- When the compressor is running, the compressed air passes through a cooler and a dryer that removes moisture from the compressed air.
- Condensate is also blown out of the drain on the dryer, which can be heard as a hiss when the compressor shuts down or when the drying chambers switch.

Risk of damage to pneumatic parts.

The pressure switch (1) was set up by the manufacturer and any other changes to its settings may be performed only by a qualified technician trained by the manufacturer.

18. SWITCHING OFF THE COMPRESSOR

- Switching off the compressor due to service or any other reason shall be performed using the pressure switch (1) by turning the switch (2) to position „0“ (Fig. 17) and **pulling out the mains plug from the socket**. This disconnects the compressor from the mains supply
- Vent the pressure in the air tank to zero by opening the drain valve.

AD DRYER

19. PRINCIPLE OF OPERATION

The dryer is controlled by a signal from the compressor pressure switch.

When the compressor is running, compressed air enters the cooler where it is cooled, and a portion of the condensed moisture is separated in the integrated cyclical separator. Solenoid valve V3 is located in the bottom of the separator, and it drains off the condensate from the separator at regular intervals.

Air then passes to the chamber of logic valve OR1, which senses the pressure. If the pressure in the valve chamber does not reach $P1 = 3$ bar, regeneration solenoid valves V1 and V2 (normally closed type) remain closed and regeneration is not initiated. If the pressure in the logic valve chamber exceeds a minimum of 3 bar, the chambers switch based on the logic described in Subsection 19.1.

When the signal from the compressor pressure switch deactivates, regeneration valves V1 and V2 open for a short time to vent the compressor aggregate and to vent the dryer chambers.

Dryer venting and regeneration is completed via the hose shared with the condensate separator drain and for this reason the condensate outlet cannot be connected directly to a drain.

19.1. Dryer regeneration

Dryer regeneration is initiated when working pressure $P1 = 3$ bar. During regeneration, valves V1 and V2 switch cyclically, whereby there is a short delay (dT) to equalise the pressures in the chambers and during which both valves are closed before each switching cycle.

Working pressure is in chamber 1 and the chamber is in the drying phase, while there is atmospheric pressure in chamber 2 and this chamber is in the regeneration phase. The chamber is regenerated by air passing through the dryer nozzle.

19.2. Energy saving mode (only for models with a dew point sensor)

Dryers equipped with a dew point sensor have additional control over the purge air so as to maintain the defined dew point.

Once the required dew point is met, regeneration is stopped to optimise the quantity of purge air. When regeneration is stopped, both valves V1 and V2 are closed (and therefore chamber purging with dry air is stopped).

They continue to alternate following their standard interval to equalise the loading on the chambers.

If the dew point increases, *Energy saving* mode is deactivated and the dryer continues to operate with regeneration activated.

Duration of deactivated regeneration mode

Energy saving mode with deactivated regeneration is time limited. One cycle with the standard regeneration of both chambers is added after the end of a defined interval (4 chamber switches).

19.3. Alarm system

All device alarms are classified as low priority technical alarms pursuant to EN 60601-1-8. Every alarm is indicated by a corresponding visual signal (see Chapter 20) and the visual signal flashes. Call an authorised service provider to conduct repairs required in the event of a failure.

Alarm states indicate a potential failure of the device. Therefore, alarm signals may not be deactivated.

All alarm signals are intermittent – when the cause of the alarm no longer exists, the alarm signals are deactivated.

Operator's position. No specific operator's position is required. The operator may remain within range of a visual warning signal depending on actual operating conditions.

Low pressure. This signal is active when the dryer pressure sensor reads low pressure in the dryer under the defined limit of 3 bar. An air leak from the compressor/dryer, malfunction of the control solenoid valves on the dryer or a software error may activate this alarm during operation.

High pressure. This signal is active when the dryer pressure sensor reads an increase in pressure over the high limit. The aggregate will then shut down and remain shut-down until the pressure decreases. If the aggregate does not shut down, excess pressure is vented by the safety valve on the compressor's air tank.

High dew point. This signal is active when the dew point sensor reads a high dew point value.


Expiration of service interval. This signal is activated until service is performed and the interval reset, or a new service interval is entered.

20. USER INTERFACE/SETTINGS

20.1. Meaning of symbols on the display

Symbol representing the user manual, or to read the user manual.


Symbol expressing the ambient temperature (from 25°C) related to the selected mode.

The program adjusts the chamber pre-filling in the given mode.

Symbol expressing the ambient temperature (from 25°C to 30°C) related to the selected mode.

The program adjusts the chamber pre-filling in the given mode.

Symbol expressing the ambient temperature (over 30°C) related to the selected mode.

The program adjusts the chamber pre-filling in the given mode.

Back to the main menu on the 4 main screens.

Go back.

Symbol for the service interval.


Symbol for when the dryer chamber is currently in regeneration.

Symbol for when the dryer chamber is currently drying.


Symbol for chamber pre-filling – equalisation of pressure in chambers.


20.2. Initial settings for the user interface

The option to select language, time and date appears when the program is first launched. The service technician selects the given language, time and date during installation based on the geographic location of the dryer installation. The program saves and remembers this setting. The option does not appear when the program is restarted. Confirming this option will also start the service hours countdown.


20.3. Setting the date and time during commissioning

It is recommended to set the date and time during commissioning. These settings are changed in the service menu in the SETUP section. Access to the service menu is protected by a code (1992). In SETUP you have to tap down arrow to set date and time. In the same way, it is possible to reset the user interface language via LANGUAGE SETUP.

20.4. Setting the drying mode during commissioning

It is recommended to set the air-drying mode during commissioning. These settings are changed in the service menu in the SERVICE SETUP section. Access to the service menu is protected by a code.

The drying mode is set based on the expected pressure range in which the device will operate and the version of the compressor (with or without a cabinet).

- compressor without a cabinet $t_0 < 25\text{ °C}$ - set to ECONOMY mode
- compressor without a cabinet $t_0 = 25\text{ °C} - 30\text{ °C}$ - set to STANDARD mode
- compressor without a cabinet $t_0 > 30\text{ °C}$ - set to ENHANCED mode
- compressor in cabinet - set to ENHANCED mode


t_0 - ambient temperature

21. MAIN SCREENS

The home screen opens once the user interface language and drying mode are set. There are 4 main screens accessed from the home screen.

The arrows at the bottom of the screen are used to navigate between the individual screens.


21.1. Home screen

(*this information is only shown if the compressor is equipped with a dew point sensor)

- Basic information and measured parameters are shown on the dryer home screen. The dryer operating hours indicator is located at the top of the screen. The home screen has 4 main tiles that display the following information:
- Dryer operating mode (see chap. 20.1)
- Current pressure in the dryer (min. working pressure to activate dryer regeneration is 3 bar)
- Current pressure dew point value (only for compressors equipped with dew point sensors)
- Remaining time to the next service interval (see chap.22.1, 22.2 – service intervals).

21.2. Dryer operation screen

The dryer operating hours indicator is located at the top of the screen. The dryer operation screen is also divided into 4 sections, which describe the dryer chamber working cycles using a graphical depiction of the operating states of each chamber.

Symbol for when the dryer chamber is currently in regeneration.

Symbol for when the dryer chamber is currently drying.

Symbol for chamber pre-filling – equalisation of pressure in chambers.

21.3. Language settings screen

LANGUAGE SETUP	
CZ	SK
GER	UK
RUS	HUN
FR	
◀	▶

The language settings screen is one of the main screens.

21.4. Service menu screen

SERVICE MENU	CODE: * * * *
SETUP	1 2
TEST	2 4
RESET	5 6
	7 8
	9 ←
◀ ▶	↩ ENTER

The service menu is intended exclusively for service personnel. Access to individual service menu items is protected by the password 1992.

22. INDICATION OF MAINTENANCE AND ALARM STATES

The adsorption dryer controller is capable of measuring selected parameters (pressure and pressure dew point in versions with a dew point sensor) and for assessing the basic dryer operating states.


Based on the measured parameters and service details, the controller displays information on faults using the following notifications and warnings.

22.1. Notification - Upcoming service interval

NOTICE !	
MODE ECO	PDP -45 °C
PRESSURE 7,3 bar	60 HRS
◀ ▶	

- 100 operating hours prior to the next service interval a **yellow notification appears** (“Notice!”). **The service tile also lights up orange.**
- **The remaining hours to the next service interval are also displayed.**
- **The warning is shown until service is completed and the service interval is reset.**


22.2. Warning - Service interval exceeded

- A red warning is shown (“Warning!”) once the service interval has lapsed. The service tile also lights up red.
- The warning is shown until service is completed and the service interval is reset.
- The hours to the next service interval are also displayed.
- A list of warnings is saved in the “Warning log”
- The dryer continues operating normally


22.3. Alarm conditions

Low/high pressure

- A red warning is shown (“Warning!”) if low pressure (<3 bar) or high pressure (> 8.5 bar, or > 10.5 bar depending on the pressure version) is detected. At the same time, the current pressure tile lights up red.
- The warning is displayed as long as the alarm condition endures.
- A list of warnings is saved in the “Warning log”
- The dryer is not operating normally, regeneration is not initiated, the compressor must be shut down.

High dew point

- This warning is only accessible on dryers with dew point sensors
- A red warning is shown (“Warning!”) if a high dew point is detected (> -20°C, or > -40°C depending on the version of the dryer). At the same time, the current pressure dew point tile lights up red.
- The warning is displayed as long as the alarm condition endures.
- A list of warnings is saved in the “Warning log”
- The dryer does not operate properly or is used outside of recommended conditions for use, the compressor must be shut down.

A high dew point alarm may appear when the unit is first placed into service. The alarm will persist until the dryer regenerates. The regeneration time of the dryer depends on many parameters of use of the device. Dryer regeneration should take no more than 24 dryer working hours. If the dryer cannot meet the dew point set-point during the first 24 hours, please contact service personnel.

22.4. “Call Service” screen

- The “Call Service” screen alternates with the main screen if a notification or warning occurs.

PRODUCT MAINTENANCE

23. PRODUCT MAINTENANCE

The operator shall ensure completion of repeated testing of the device at least once every 24 months (EN 62353) or at intervals defined by applicable national legal regulations. A record of these test results shall be completed (e.g. per EN 62353, Annex G) together with the measurement methods.

The device has been designed and manufactured to keep its maintenance to a minimum. For correct and reliable operation of the compressor perform the following operations:

Risk of servicing by persons without required qualification.

Repair works beyond standard maintenance (see chapter 23.1) may only be performed by a qualified technician (organization authorized by the manufacturer) or manufacturer's customer service.

Standard maintenance works (see chapter 23.1) may only be carried out by a trained operator.

Only use spare parts and accessories approved by the manufacturer.

Risk of injury and damage to the device.

Before any maintenance work, it is necessary to:

- check, if it is possible to disconnect compressor from the supplied equipment, so that there is no risk of damage to a person using the given appliance or any other material damages
- switch off the compressor
- disconnect it from the mains (pull out the mains plug)
- vent the compressed air from the air tank

Risk of injury when venting the compressed air.

Protect eyesight – wear goggles when venting the compressed air from the compressor pneumatic system (air tank).

Risk of burn.

Pump components (head, cylinder, and pressure hose) have high temperature during and shortly after compressor operation – do not touch these components!

Let the device cool before any product maintenance, service or connection/disconnection of pressurized air!

The removed grounding conductor during service must be connected back to the original position after completing the service.

- For any maintenance or repair work, the compressor may be removed from the cabinet by pulling it from the cabinet (on rolling casters as much as the cord between the compressor and the cabinet allows). See also chap. 23). The required maintenance or repair work may then be performed.
- If the cord between the compressor and the cabinet becomes disconnected (by pulling the mains plug, chap. 13.2), the compressor shall not be functional. In order to bring the compressor into function again, please follow chap. 13.2 or 23.10 or 23.11.

23.1. Maintenance intervals

Time interval ^{b)}	50 Hz		60 Hz		Once a day	Once a week	Once a year	Once every 2 years	4000 hours	3200 hours	4800 hours	6000 hours	8000 hours	10000 hours	12000 hours	16000 hours	20000 hours	24000 hours	Chap.	Spare parts set	Performed by		
																					operator	qualified technician	
Check of product operation	x																		23.2	-		operator	
Check of product function						x													9	-		operator	
Check of pneumatic connections leakage and device inspection							x												23.3	-		qualified technician	
Check of electrical connection							x												23.4	-		qualified technician	
Check of fan and cooler							x												23.8	-		qualified technician	
Check of safety valve							x												23.6	-		qualified technician	
Replacement of pump inlet filter and pre-filter ^{a)}							x		x				x						23.7	025200139-000 025200150-000		qualified technician	

^{a)} data is in hours, if not available, data is in years

^{b)} time interval reduces by 20 % for compressor variants 60 Hz

(2000h./1600h., 4000h./3200h., 6000h./4800h., 8000h./6400h., 10000h./8000h., 12000h./9600h.)

^{c)} Whichever occurs first: 4 years/8,000 hours.

Time interval ^{b)}	50 Hz		60 Hz		Once a day	Once a week	Once a year	Once every 2 years	4000 hours	6000 hours	8000 hours	10000 hours	12000 hours	16000 hours	20000 hours	24000 hours	Chap.	Spare parts set	Performed by
	3200 hours	4800 hours	6400 hours	8000 hours															
Replacement of cassettes with adsorbent media ^{c)}																	See service manual	603022476-000	qualified technician
Replacement of the dryer's logic valve ball									x								See service manual	069000317-000	qualified technician
Replacement of the dryer's solenoid valve ^{d)}																	See service manual	025300117-000	qualified technician
Dew point sensor calibration																	23.9	-	qualified technician
Conduct a "repeated test" per EN 62353																	23	-	qualified technician

a) data is in hours, if not available, data is in years

b) time interval reduces by 20 % for compressor variants 60 Hz

(2000h./1600h., 4000h./3200h., 6000h./4800h., 8000h./6400h., 10000h./8000h., 12000h./9600h.)

c) Whichever occurs first: 4 years/8,000 hours

d) Whichever occurs first: 8 years/20,000 hours

23.2. Check of product operation

- Check pump condition – the pumps must operate evenly without excessive vibrations or noise. If there is any negative result, find the cause of the given condition or call service engineers.
- Check fans operation (visually) – the fans must be operating when the fans are running. If there is a negative result, look for the cause of the given condition or call service engineers.
- Check if the power cord and pneumatic hoses are not damaged. Replace the damaged parts or call service engineers.
- Check ambient temperature – ambient temperature must be below temperature limit (40°C). Improve cooling in the room if the temperature is higher.
- Open the plug on the condensate collecting bottle and drain the condensate.

23.3. Check of pneumatic connections leakage and device inspection

Leakage check

- Check pneumatic connections of the compressor when the compressor is operating – pressurizing.
- Use a leakage analyser or soapy water to check all connections for leaks. Tighten or release the connection where leaks are found.

Device inspection

- Check condition of the compressor pump – even operation, appropriate noise level.
- Check fans operation – fans must be running in the defined cycles of the compressor operation.
- Check function of the temperature switch (B2) – heat the temperature switch to temperature above 40°C (e.g. using a heat gun – Be careful, do not operate with high temperature on plastic components in the surrounding, it can cause deformation of plastic materials). Fan E1 starts up once the temperature reaches 40°C – the compressor must be under voltage.
- Check filters condition – filters must be undamaged and sufficiently clean.
- Check condition of the pump itself, check if there is no dirt inside the crankcase or clearance in the crankshaft.
- Replace any defective parts if there are any failures.

23.4. Check of electrical connections

Risk of electric shock.

Check electrical connections on the device disconnected from the mains.

- Check mechanical function of the main switch.
- Check if the power cable and conductors are not damaged.
- Visually check if cables are connected to the terminal box.
- Check all screw connections of the green-yellow PE grounding conductor.
- Inspection of electrical connections on the AD dryer – see the service manual.

23.5. Condensate drain

A wet floor resulting from overflow from the vessel poses a slip hazard.

Regularly empty the condensate collection tank.

Injury hazard.

Do not connect the condensate drain directly to a waste drain!

Passers-by may be injured.

Switch off the compressor every time before emptying the tank!

Injury hazard.

Carefully and slowly open the drain valve. Quickly opening the drain valve is accompanied by a loud sound and an uncontrollable ejection of the accumulated condensate.

Condensate is automatically drained into the condensate collection vessel.

- Monitor the level in the vessel using the 1 l or 2 l markings (depending on the volume of the vessel), and empty at least once a day

Fig. 18: Check of condensate collection vessel

If necessary, it is possible to drain the air tank using the drain valve (1) from the lower side of the air tank (Fig. 19).

Switch off the compressor at the mains and reduce the air pressure in the device to a maximum of 1 bar, e.g. by venting air through the connected equipment.

Use a hose to drain with the free end placed in a suitable container and then open the valve to drain condensate from the tank.


Fig. 19: Condensate drain

Before the following checks it is required:

- **Compressor with cabinet** – release the door lock and open the cabinet door.

23.6. Check of safety valve

- Turn the screw (2) on the safety valve (1) several times to the left until the safety valve releases the air.
- Let the safety valve vent for a few seconds.
- Turn the screw (2) fully to the right, the valve must be closed now.


Fig. 20: Check of safety valve

Risk of pressure increase when safety valve is damaged.

Safety valve must not be used to depressurize the air tank. It could damage the safety valve function. The valve is set to the maximum allowed pressure, it is tested and marked.

Safety valve setting must not be adjusted.

Risk of injury when compressed air is vented.

Protect your eyesight when checking the safety valve – wear protective goggles.

23.7. Inlet filter and pre-filter replacement

The inlet filter (1) is located inside the compressor crankcase cover.

Inlet filter replacement:

- Pull out the rubber plug (2) by hand.
- Remove the dirty inlet filter (1).
- Insert a new filter and mount the rubber plug.

Pre-filter replacement:

- Pull out the pre-filter by hand (3).
- Replace the element and re-insert.


Fig. 21: Inlet filter and pre-filter replacement

23.8. Check of fan and cooler

The entire device and especially the compressor fan, cooler fan and cooler must be kept clean (Fig. 1) – vacuum or blown down the cooling fins and fans with compressed air to remove any dust from the surface.

23.9. Dew point sensor calibration

It is recommended to calibrate the dew point sensor every 2 years on compressors equipped with this sensor. Calibration is performed by the manufacturer (Vaisala) and must be ordered online.

Sensor calibration ordering process:

- Visit the website www.vaisala.com, and then navigate to the *Services* section and select *Calibration and repair services*.
- Then enter the type of sensor (sensor type: **DMT143**) and select a calibration type, then follow the provided instructions.

23.10. Procedure for connecting a compressor disconnected from the cabinet

Prior to any maintenance or repair work, switch off the compressor and disconnect it from the mains (pull out the mains plug).

The compressor (without cabinet) requires that the jumper is always mounted in the terminal strip for correct operation (Fig. 22, Pos. B). It then replaces the switch on the cabinet.

If there is no jumper in the terminal strip, the compressor will not work!

If the cord between the compressor and the cabinet becomes disconnected (by pulling the mains plug) and the compressor is removed from the cabinet, the compressor will not be functional. Therefore, it is necessary to first create a connection to the terminal strip with a jumper (replacing the circuit breaker function) (Fig. 22).


PROCESS:

Jumper installation (A-B process):

- Disconnect the product from the mains by pulling the plug out of the socket.
- Remove the cover on the electrical panel (on the compressor).
- The jumper is not in the terminal strip – A.
- Insert the jumper into the terminal strip – B.
- Reinstall the cover on the electrical panel.
- Connect the device to the electrical mains.
- Activate the compressor by turning on the switch on the pressure switch.

Compressor - when operating outside the cabinet, the electrical connection with an installed jumper **must** be created (Fig. 22, Pos. B).

Fig. 22

23.11. Procedure for connecting a compressor to a new cabinet

Prior to any maintenance or repair work, switch off the compressor and disconnect it from the mains (pull out the mains plug).

The compressor in a cabinet requires that the jumper is not mounted in the terminal strip for correct operation (Fig. 22, Pos. A). The switch on the cabinet switches the entire product on and off, including the compressor.

If there is a jumper in the terminal strip, the switch on the compressor cabinet will not work!

To connect a new cabinet to a compressor that previously worked without a cabinet, it is necessary

to cancel the connection in the compressor terminal strip by performing the procedure as follows. (Fig. 22). See also chap. 13.1 and 13.2.

PROCESS:

Jumper removal (B-A process):

- Disconnect the product from the mains by pulling the plug out of the socket.
- Remove the cover on the electrical panel (on the compressor).
- The jumper is in the terminal strip – B
- Remove the jumper from the terminal strip – A
- Reinstall the cover on the electrical panel.
- Insert the compressor into the cabinet.
- Connect the cabinet to the compressor by inserting the power cord with connector (1) into a power socket (2). (Fig. 11)
- Disconnect the display on the AD dryer and connect the display at the front door to the cabinet (see the service manual).
- Connect the device to the electrical mains
- Activate the compressor by turning on the switch on (2) the pressure switch (1) and the switch (5) on the cabinet (Fig. 17)

Compressor with cabinet - when operating in the cabinet, there **must not** be an electrical connection. (Fig. 22, Pos. A)

Note: **In case of using a compressor without a cabinet or in case of maintenance, a jumper which is not installed into the terminal strip should be kept** in a designated place next to the electrical panel.

23.12. Cleaning of the exterior surfaces of the product

Use neutral agents for cleaning of the external surfaces of the product.

The use of aggressive cleaning agents and comprising alcohol and chlorides may lead to the damage of the surface and the discoloration of the product.

24. LONG-TERM SHUTDOWN

If the compressor is not going to be used for long period, it is recommended to drain all condensate from the air tank and switch on the compressor for about 10 minutes, keeping the drain valve open. Thereafter switch off the compressor using the switch (2) on the pressure switch (1) (Fig. 17), close the drain valve and disconnect the device from the mains.

25. DISPOSAL OF DEVICE

- Disconnect the device from the mains.
- Release air pressure in the pressure tank by opening the drain valve.
- Dispose of the device according to the applicable local rules.
- Order a specialized organization to sort and dispose of waste.
- Product components after its operational lifetime have no negative environmental effect.

TROUBLESHOOTING**Risk of electric shock.****Before any of the following operations on the device, disconnect the device from the mains (pull out the mains plug).****Risk of injury during work with pneumatic components under pressure.****Before any of the following operations on the device, it is necessary to decrease pressure in the air tank and in the pneumatic system to zero.****Any operations concerning the troubleshooting may be performed only by a qualified technician.**

Malfunction	Possible cause	Solution
Compressor does not switch on	No voltage in the pressure switch	Check voltage in socket
		Check circuit breaker switch – switch to position switched-on „I“
		Loose conductor from terminal - repair
		Check electrical cord – replace defective cord
	Motor winding failure, damaged thermal protection	Replace motor or windings
	Capacitor failure	Replace capacitor
	Seized piston or other rotating part	Replace damaged components
	Pressure switch does not switch	Check function of pressure switch
The cabinet switch does not switch off the compressor	There is a compressor in the cabinet that previously worked separately (there is a jumper in the compressor terminal strip)	Remove the jumper from the terminal strip, see chap. 23.11
Compressor often switches on	Air leakage in pneumatic system	Check pneumatic system – seal loose connections
	Non-return valve (SV) leakage	Clean non-return valve, replace seals, replace non-return valve
	Large volume of condensed liquid in pressure tank	Drain condensed liquid
	Low compressor capacity	Check time of filling the air tank
Low pressure in the air tank (compressor running constantly)	High air consumption of supplied equipment	Decrease air consumption Use compressor with higher capacity
	Leakage in pneumatic system	Check pneumatic system – seal loose connections
	Low capacity of the pump	Clean / replace the pump
	Pump failure	Clean / replace the pump
	Dryer failure	Replace dryer
Prolonged operation of the compressor	Air leakage in pneumatic system	Check pneumatic system – seal loose connection
	Worn piston ring	Replace worn piston ring
	Dirty inlet filter	Replace dirty filter with new filter

	Incorrect function of solenoid valve	Repair or replace fan or coil
Compressor is noisy (knocking, metal noises)	Damaged piston bearing, piston rods, motor bearing	Replace damaged bearing
	Loose (cracked) dampening element (spring)	Replace damaged spring
Dryer does not dry (condensed water in the air) *	Defective cooler fan	Replace fan
	Damaged dryer	Check power source Replace dryer
	Defective automatic condensate drain	Clean / replace
Degraded drying performance – high-pressure dew point (condensed water in the air)	Low operating pressure	Reduce the demand for air, check the output from the compressed air source, fix any leaks in the distribution system
	Regeneration solenoid valve not working	Check coil operation, replace if damaged Inspect the condition of the valve – clean the valve or replace if problems persist
	Air regeneration nozzle plugged	Clean or replace the nozzle (see product maintenance)
	Cooling fan not working	Check the power source to the fan Replace damaged fan
	Dirty cooler	Inspect the cooler and clean as necessary
	Silencer plugged at outlet from regeneration valve	Inspect the silencers. Clean or replace the silencer if flow resistance is too high or if heavily soiled.
Dryer emitting high levels of noise	Damaged fan	Replace damaged fan
	Damaged silencer	Replace the silencer
	Air leaking through relief valve at dryer inlet	Check the dryer connection to the mains and dryer connections, check the dryer operation, check the dryer's working pressure, and replace defective components.
Air leaking through relief valve at dryer inlet	High working pressure from air source	Check the pressure setting on the compressed air source
	Dryer inlet solenoid valve not working	Check the coil operation, replace if damaged Inspect the condition of the valve – clean the valve or replace if problems persist
	High pressure in equipment resulting from plugged filters	Check the internal filters and accessory filter assemblies. Clean or replace dirty filters.

)* Clean inside surfaces of the air tank thoroughly and remove all condensed liquid after dryer failure.

Once a dryer fault is cleared and after its reassembly, the dryer should quickly be regenerated, best when using continuous compressor operation at a pressure of around 1 bar for a period of at least 1 hour; then check the dryness of the compressed air.

To protect the connected equipment from any damage, it is necessary to check humidity of the released air from the air tank (see chapter Technical data)

26. REPAIR SERVICE

Guarantee and post-guarantee repairs are provided by the manufacturer or organizations and technicians approved by the manufacturer.

Warning

The manufacturer reserves the right to make changes on the device, which will not significantly affect properties of the device.

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1. Product: (model) DK50 2V/M DK50 2V S/M		2. Serial number:	
3.1. User's name:			
3.2. Address of installation:			
4. Equipment connected to the compressor:			
5. Installation / Commissioning:		6. Contents of operator training:	
Product completeness check **	Y	Description of the product and functions**	Y
	N		N
Documentation completeness check **	Y	Product operation: turning on/off, controls, control procedures, data on the display panel, alarms, operation in alarm conditions**	Y
	N		N
Installation/connection to equipment **	Y	Product maintenance: maintenance intervals, maintenance procedure, service intervals, operating activities**	Y
	N		N
Functional test **	Y	Safety measures, warnings – their meaning and compliance **	Y
	N		N
Notes:			
7. Operator instructed on safety measures, operations and maintenance:			
Name:		Signature:	
Name:		Signature:	
Name:		Signature:	
8. Installation and instruction performed by: First name/Last name		Signature:	
Company:		Address:	
Phone:			
Email:		Date:	
9. Distributor:			
Company:		Address:	
Contact person:			
Phone:		Email:	

** mark with an "X" in points 5 and 6 (Y - yes /N - no). Enter any observations from points 5 and 6 into the "Notes" section

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