

DK50 PLUS/M

DK50 2V/M

- EN USER MANUAL
- DE BENUTZERHANDBUCH
- FR MANUEL DE L'UTILISATEUR
- RU РУКОВОДСТВО ПОЛЬЗОВАТЕЛЯ
- PL INSTRUKCJA OBSŁUGI
- SK NÁVOD NA POUŽITIE
- CS NÁVOD K POUŽITÍ

CE 2460

COMPRESSOR
KOMPRESSOR
COMPRESSEUR
КОМПРЕССОР
SPRĘŻARKA
KOMPRESOR

DK50 PLUS/M
DK50 2V/M

EKOM spol. s r. o.
Priemyselná 5031/18
SK-921 01 Piešťany
Slovak Republic
tel.: +421 33 7967255
fax: +421 33 7967223

www.ekom.sk
email: ekom@ekom.sk

DATE OF LAST REVISION
DATUM DER LETZTEN ÜBERARBEITUNG
DATE DE LA DERNIÈRE RÉVISION
ДАТА ПОСЛЕДНЕЙ РЕДАКЦИИ
DATA OSTATNIEJ AKTUALIZACJI
DÁTUM POSLEDNEJ REVÍZIE
DATUM POSLEDNÍ REVIZE

01/2025

NP-DK50 PLUS, 2V-ADS-
6_01-2025
112000583-000

EN	<u>CONTENTS</u>	<u>5</u>
DE	<u>INHALT</u>	<u>51</u>
FR	<u>TABLES DES MATIÈRES</u>	<u>99</u>
RU	<u>ОГЛАВЛЕНИЕ</u>	<u>147</u>
PL	<u>SPIIS TREŚCI</u>	<u>199</u>
SK	<u>OBSAH</u>	<u>247</u>
CS	<u>OBSAH</u>	<u>293</u>

CONTENTS

IMPORTANT INFORMATION	6
1. CONFORMITY WITH THE REQUIREMENTS OF EUROPEAN UNION DIRECTIVES	6
2. INTENDED USE	6
3. CONTRAINDICATIONS AND SIDE-EFFECTS	6
4. SYMBOLS.....	6
5. WARNINGS	7
6. STORAGE AND TRANSPORT CONDITIONS	9
PRODUCT DESCRIPTION	10
7. VARIANTS	10
8. ACCESSORIES	10
9. PRODUCT FUNCTION	11
TECHNICAL DATA	15
INSTALLATION	22
10. INSTALLATION CONDITIONS	22
11. COMPRESSOR ASSEMBLY	22
12. PNEUMATIC CONNECTION	24
13. ELECTRICAL CONNECTION	26
14. PLACEMENT OF THE COMPRESSOR IN THE CABINET	28
15. COMMISSIONING	30
16. PNEUMATIC AND ELECTRICAL DIAGRAMS	31
OPERATION	34
17. SWITCHING ON THE COMPRESSOR.....	35
18. SWITCHING OFF THE COMPRESSOR	35
ADS DRYER	36
19. PRINCIPLE OF OPERATION	36
20. BASIC SCREENS	36
21. SERVICE SCREENS	38
PRODUCT MAINTENANCE	40
22. PRODUCT MAINTENANCE.....	40
23. LONG-TERM SHUTDOWN.....	47
24. DISPOSAL OF DEVICE	47
TROUBLESHOOTING	48
25. REPAIR SERVICE	50

IMPORTANT INFORMATION

1. CONFORMITY WITH THE REQUIREMENTS OF EUROPEAN UNION DIRECTIVES

This product conforms to the requirements of the Medical Device Directive (93/42/EEC), and is safe for the intended use if all safety instructions are followed.

2. INTENDED USE

The compressor is used as a source of clean, oil-free compressed air to power dental unit, instruments, and devices in laboratories where the parameters and properties of the compressed air are suitable for the specific application.

Compressed air supplied by the compressor is unsuitable for use with artificial lung ventilation devices.

Any other use of the product beyond its intended use is considered an incorrect use. The manufacturer is not liable for any damages or injuries resulting from the incorrect use.

3. CONTRAINDICATIONS AND SIDE-EFFECTS

There are no contraindications or side-effects known.

4. SYMBOLS

The following symbols and marks are used in the User manual, on the device and its packaging:

General warning

Warning - risk of electric shock

Warning - compressor is controlled automatically

Warning - hot surface

General caution

Refer to instruction manual

CE – marking

Medical device

Serial number

Article number

Unique Device Identifier

	Swiss Authorised Representative
	Swiss Importer
	Protecting earthing
	Terminal for ground connection
	Fuse
	Package handling label – FRAGILE
	Package handling label – THIS SIDE UP
	Package handling label – KEEP DRY
	Package handling label – TEMPERATURE LIMITS
	Package handling label – LIMITED STACKING
	Package label – RECYCLABLE MATERIAL
	Manufacturer

5. WARNINGS

The product is designed and manufactured to be safe for the user and the surrounding environment when used in the defined manner. Keep the following warnings in mind.

5.1. General warnings

READ THE USER MANUAL CAREFULLY BEFORE USING THE DEVICE AND KEEP IT FOR FUTURE USE!

- The user manual supports the correct installation, operation and maintenance of the product. Careful review of this manual will provide the information necessary for the correct operation of the product for its intended use.
- Original packaging should be retained for the return of the device. Only the original packaging ensures protection of the device during transport. If it is necessary to return the product during the warranty period, the manufacturer is not liable for damages caused by improper packaging.
- This warranty does not cover damages originating from the use of accessories other than those specified or recommended by the manufacturer.
- The manufacturer only guarantees the safety, reliability and function of the device if:
 - installation, new settings, changes, expansion, and repairs are performed by the manufacturer or an organisation authorised by the manufacturer

- the product is used pursuant to the user manual
- The user manual corresponds to the configuration of the product and its compliance with the applicable safety and technical standards at the time of printing. The manufacturer reserves all rights for the protection of its configuration, methods and names.
- Translation of the user manual is performed in accordance with the best available knowledge. The Slovak version is to be used in the event of any uncertainties.
- This user manual is the original instructions. Translation is performed in accordance with the best available knowledge.

5.2. General safety warnings

The manufacturer designed and manufactured the product to minimise all risks when used correctly for the intended use. The manufacturer considers it its obligation to lay down the following general safety precautions.

- Use and operation of the product must comply with all laws and local regulations valid in the place of use. The operator and user are responsible for following all the appropriate regulations in the interests of performing work safely.
- Only the use of original parts guarantees the safety of operating personnel and the flawless operation of the product itself. Only the accessories and parts mentioned in the technical documentation or expressly approved by the manufacturer should be used.
- The operator must ensure that the device is functioning correctly and safely before every use.
- The user must be familiar with the operation of the device.
- Do not use the product in environments with a risk of explosion.
- The user must inform the supplier immediately if any problem directly related to the operation of the device occurs.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

5.3. Safety warnings on protection from electric current

- The device must only be connected to a properly installed, earthed socket.
- Before the product is plugged in, ensure that the mains voltage and frequency stated on the product are in compliance with the values of the mains.
- Prior to putting the device into operation check for any damage to the connected pneumatic lines and electrical wiring. Replace damaged pneumatic lines and electrical wirings immediately.
- Immediately disconnect the product from the mains (remove power cord from the socket) in hazardous situations or in the case of a technical malfunction.
- During all repairs and maintenance, ensure that:
 - the mains plug is removed from the power socket
 - pressure is vented from the air tank and pipes
- The product shall only be installed by a qualified technician.

6. STORAGE AND TRANSPORT CONDITIONS

The compressor is shipped from the manufacturer in transport packaging. This protects the product from damage during transport.

Potential for damage to pneumatic components.

The compressor may only be transported when all air pressure has been vented. Before moving or transporting the compressor, release all air pressure from the tank and hoses and drain condensate from the air tank.

Keep the original factory packaging in case the device needs to be returned Use the original factory packaging during transport as it provides optimum protection for the product. If it is necessary to return the product during the warranty period, the manufacturer is not liable for damages caused by improper packaging.

The compressor is shipped in a vertical position and must be secured using transport straps.

Protect the compressor from humid and dirty environments and extreme temperatures during transport and storage. Do not store near any volatile chemical substances.

If not, please dispose of the original packaging material in an environmentally-friendly way. The packaging cardboard can be recycled with old paper.

Storing or shipping the equipment in any conditions other than those specified below is prohibited.

6.1. Ambient conditions

Products may only be stored and transported in vehicles that are free of any traces of volatile chemicals under the following climactic conditions:

Temperature

–25°C to +55°C

Relative humidity

max. 90% (non-condensing)

PRODUCT DESCRIPTION

7. VARIANTS


The compressor is manufactured according to its intended application in the following variants:

DK50 PLUS/M Compressor on the base with adsorption air dryer


DK50 PLUS S/M Compressor in cabinet with adsorption air dryer

DK50 2V/M Compressor on the base with adsorption air dryer

DK50 2VS/M Compressor in cabinet with adsorption air dryer

DK50 PLUS/M

DK50 2V/M

DK50 PLUS S/M
DK50 2V S/M

8. ACCESSORIES

Accessories that are not included in the standard order must be ordered separately.

Set of compressed air outlet filters

The compressor may be equipped with a set of filters if specified. The filter set may be equipped with an air pressure regulator.

Where a different level of air filtration is required, this requirement must be agreed upon with the supplier and specified in the order.

Type	Use	Level of filtration / μm /	Article number
FS 30F		1	604014119-014
FS 30M	DK50 PLUS/M	1+0.1	604014119-015
FS 30S	DK50 2V/M	1+0.01	604014119-017
FS 30A		1+0.01+A**	604014119-018

**) A –submicrofilter with activated carbon

Filter set regulator assembly

The compressor may be equipped with a pressure regulator of the compressed air outlet if specified. The regulator must be selected according to the application to the filter set, or separately. The regulator shall ensure constant pressure at the outlet.

Type	Use	Article number
Regulator complete (to the filter assembly)	DK50 PLUS/M	604014130-000
REG11 (separate)	DK50 2V/M	447000001-077

Filter set brackets

A suitable bracket must be ordered for every filter set.

Type	Use	Article number
Compressor-mounted bracket	DK50 PLUS/M DK50 2V/M	603014134-000

Equipotential bonding socket

The socket allows an equipotential bonding

Type	Use	Name	Article number
POAG-KBT6-EC	DK50 PLUS S/M	Mains plug connector	033400075-000
FLEXI-S/POAG-HK6	DK50 2VS/M	Conductor (1 m)	034110083-000

9. PRODUCT FUNCTION

9.1. Compressor with adsorption dryer

The air pump (1) draws in air through an inlet filter (8) and compresses it, feeding it to the cooler (10) through the check valve (3) into the dryer (9), and feeds the clean and dry air into the air tank (2). The dryer continuously dries the compressed air. Condensate is drained from the air tank via the drain valve (7), depending on the drying controls. Compressed, oil-free filtered air is stored in the air tank ready for further use.

It is not necessary to drain the air tank.

A portion of the dried air from the air tank is then used to regenerate the adsorption dryer during a break in the compressor operation.

Risk of damage to the device.

A compressor with a type ADS adsorptive dryer may only be operated in the S3 operating mode. The compressor is optimised for the S3-70% operating mode. It is not recommended to use the compressor in a mode higher than S3-80%.

9.2. Compressor cabinet

The soundproof cabinet provides compact covering of the compressor, dampens noise efficiently, providing sufficient exchange of cooling air. With its design It can be placed in a dentist's office as a part of furniture. The fan (10) under the compressor pump provides cooling of the compressor. It is in operation at the same time with the compressor motor or after switching on the temperature switch at temperature above 40°C. After cooling the cabinet area approximately under 32°C, the fan turns off automatically.

Door of the casing with right opening may be changed to left opening (see chap. 14.1).

Risk of compressor overheating.

Make sure that there are no obstacles at the cooling air inlet into the cabinet (around the bottom part of the cabinet) and at the hot air outlet on the top back side of the cabinet.

If the compressor is placed on a soft floor, e.g. carpet, create space between the base and the floor or the cabinet and the floor, e.g. underlay the footings with hard pads to ensure sufficient cooling of the compressor.

Descriptions for figures 1 - 2

1 Air pump	16 Cooler fan
2 Air tank	17 Cabinet fan
3 Non-return valve	18 Rectifying screw
4 Pressure switch	19 Lock
5 Safety valve	20 Door pin
6 Pressure gauge	21 Wall stop
7 Drain valve	22 Door hinge
8 Inlet filter	23 Switch
9 Adsorption dryer	24 Cabinet pressure gauge
10 Cooler	25 Handle
11 Compressor fan	26 Supporting bracket
12 Magnetic holder	27 Casters
13 Condensate collection vessel	28 Pressure gauge hose
14 Compressed air outlet	29 Power cord
15 Circuit breaker switch	30 Dryer control electronics

Fig. 1: DK50 2V/M – Compressor with adsorption dryer



Fig. 2: Cabinet

TECHNICAL DATA


Compressors are designed for operation in dry, ventilated and dust-free indoor rooms under the following climactic conditions:

Temperature	+5°C to +40°C
Relative humidity	max. 70%

Working pressure 6 – 8 bar		DK50 PLUS/M	DK50 PLUS S/M
Nominal voltage, Frequency ^{a)}	V, Hz	230, 50/60	230, 50/60
Capacity at 6 bar (FAD)	l/min	70 ^{d)} / 80 ^{e)}	70 ^{d)} / 80 ^{e)}
Working pressure ^{b)}	bar	6.0 – 8.0	6.0 – 8.0
Rated current	A	4.2/5.2	4.3/5.3
Motor power	kW	0.55	0.55
Air tank volume	l	25	25
Air quality - filtration	µm	-	-
Maximum operating pressure of safety valve	bar	11.0	11.0
Noise level at 5 bar (LpA)	dB	≤ 63/65	≤ 50/51
Operating mode	%	S3-70%	S3-70%
PDP drying performance at 7 bar	°C	≤ -20	≤ -20
Time to fill air tank from 0 to 7 bar	s	160/140	160/140
Dimensions (net) w x d x h	mm	460x565x710	560x690x875
Net weight ^{c)}	kg	56	95
Classification under EN 60601-1	Class I.		
Classification pursuant to MDD 93/42 EEC,2007/47 EC	IIa		

- a) Specify the compressor version when ordering
- b) Consult any other range of pressure with the supplier
- c) Weight is indicative and only applies to the product without accessories
- d) In order to maintain the declared dewpoint, the maximum continuous air consumption of the connected device shall not exceed 40 l/min.
- e) In order to maintain the declared dewpoint, the maximum continuous air consumption of the connected device shall not exceed 46 l/min.


Dependence of compressor output on working pressure

Working pressure 8 – 10 bar		DK50 PLUS/M	DK50 PLUS S/M
Nominal voltage, Frequency ^{a)}	V, Hz	230, 50/60	230, 50/60
Capacity at 8 bar (FAD)	l/min	60 ^{d)} / 70 ^{e)}	60 ^{d)} / 70 ^{e)}
Working pressure ^{b)}	bar	8.0 – 10.0	8.0 – 10.0
Rated current	A	4.3/5.3	4.4/5.4
Motor power	kW	0.55	0.55
Air tank volume	l	25	25
Air quality - filtration	µm	-	-
Maximum operating pressure of safety valve	bar	11.0	11.0
Noise level at 5 bar (LpA)	dB	≤ 63/65	≤ 50/51
Operating mode	%	S3-70%	S3-70%
PDP drying performance at 7 bar	°C	≤ -20	≤ -20
Time to fill air tank from 0 to 7 bar	s	160/140	160/140
Dimensions (net) w x d x h	mm	460x565x710	560x690x875
Net weight ^{c)}	kg	56	95
Classification under EN 60601-1	Class I.		
Classification pursuant to MDD 93/42 EEC, 2007/47 EC	IIa		

- a) Specify the compressor version when ordering
- b) Consult any other range of pressure with the supplier
- c) Weight is indicative and only applies to the product without accessories
- d) In order to maintain the declared dewpoint, the maximum continuous air consumption of the connected device shall not exceed 35 l/min.
- e) In order to maintain the declared dewpoint, the maximum continuous air consumption of the connected device shall not exceed 41 l/min.


Dependence of compressor output on working pressure

Working pressure 6 – 8 bar		DK50 2V/M	DK50 2V S/M
Nominal voltage, Frequency ^{a)}	V, Hz	230, 50	230, 50
Capacity at 6 bar (FAD)	l/min	135 ^{d)}	135 ^{d)}
Working pressure ^{b)}	bar	6.0 – 8.0	6.0 – 8.0
Rated current	A	8.7	8.9
Motor power	kW	1.2	1.2
Air tank volume	l	25	25
Air quality - filtration	µm	-	-
Maximum operating pressure of safety valve	bar	11.0	11.0
Noise level at 5 bar (LpA)	dB	≤ 69	≤ 55
Operating mode	%	S3-70%	S3-70%
PDP drying performance at 7 bar	°C	≤ -20	≤ -20
Time to fill air tank from 0 to 7 bar	s	80	80
Dimensions (net) w x d x h	mm	460x565x710	560x690x875
Net weight ^{c)}	kg	61	99
Classification under EN 60601-1	Class I.		
Classification pursuant to MDD 93/42 EEC,2007/47 EC	IIa		


- a) Specify the compressor version when ordering
- b) Consult any other range of pressure with the supplier
- c) Weight is indicative and only applies to the product without accessories
- d) In order to maintain the declared dewpoint, the maximum continuous air consumption of the connected device shall not exceed 78 l/min.

Dependence of compressor output on working pressure

Working pressure 8 – 10 bar		DK50 2V/M	DK50 2V S/M
Nominal voltage, Frequency ^{a)}	V, Hz	230, 50	230, 50
Capacity at 8 bar (FAD)	l/min	110 ^{d)}	110 ^{d)}
Working pressure ^{b)}	bar	8.0 – 10.0	8.0 – 10.0
Rated current	A	8.5	8.7
Motor power	kW	1.2	1.2
Air tank volume	l	25	25
Air quality - filtration	µm	-	-
Maximum operating pressure of safety valve	bar	11.0	11.0
Noise level at 5 bar (LpA)	dB	≤ 69	≤ 55
Operating mode	%	S3-70%	S3-70%
PDP drying performance at 7 bar	°C	≤ -20	≤ -20
Time to fill air tank from 0 to 7 bar	s	80	80
Dimensions (net) w x d x h	mm	460x565x710	560x690x875
Net weight ^{c)}	kg	61	99
Classification under EN 60601-1	Class I.		
Classification pursuant to MDD 93/42 EEC,2007/47 EC	IIa		

- a) Specify the compressor version when ordering
- b) Consult any other range of pressure with the supplier
- c) Weight is indicative and only applies to the product without accessories
- d) In order to maintain the declared dewpoint, the maximum continuous air consumption of the connected device shall not exceed 64 l/min.

Dependence of compressor output on working pressure

FAD correction of capacity for altitude

Capacity given in the form of FAD („Free Air Delivery“) applies to the following conditions:

Altitude	0 m.n.m.	Temperature	20°C
Atmospheric pressure	101325 Pa	Relative humidity	0%

To calculate FAD compressor capacity in dependence on altitude, it is necessary to apply correction factor according to the following table:

Altitude [m.n.m.]	0 -1500	1501 - 2500	2501 - 3500	3501 - 4500
FAD correction factor	1	0.80	0.71	0.60

Electromagnetic compatibility declaration

Medical equipment needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided below.

Guidance and manufacturer's declaration - electromagnetic emissions		
Pursuant to IEC 60601-1-2:2014/AMD1:2020 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances		
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The equipment uses RF energy only for its internal functions. Therefore, the RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The equipment is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	The equipment is not likely to cause any flicker, as the current flow is approx. constant after the start up.	

Guidance and manufacturer's declaration - electromagnetic immunity

Pursuant to IEC 60601-1-2:2014/AMD1:2020 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances

The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.

Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV 100 kHz repetition frequency Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV L-N ±2 kV L-PE; N-PE Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruption, and voltage variations on power supply input lines IEC 60601-4-11	$U_T=0\%$, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°) $U_T=0\%$, 1 cycle $U_T=70\%$ 25/30 cycles (at 0°) $U_T=0\%$, 250/300 cycles	$U_T=>95\%$, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°) $U_T=>95\%$, 1 cycle $U_T=70\%$ (30% dip U_T), 25(50Hz)/30(60Hz) cycles (at 0°) $U_T=>95\%$, 250(50Hz)/300(60Hz) cycle	Mains power quality should be that of a typical commercial or hospital environment. The device stops and restarts automatically at each dip. In this case does not occur unacceptable pressure drop.
Power frequency (50/60 Hz) IEC 61000-4-8	N/A	N/A	The equipment does not contain magnetically sensitive components or circuits.
Radiated fields in close proximity per IEC 61000-4-39	N/A	N/A	The equipment does not contain magnetically sensitive components or circuits.

NOTE: U_T is the A.C. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity			
Pursuant to IEC 60601-1-2:2014/AMD1:2020 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances			
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d=1,2\sqrt{P}$ $d=1,2\sqrt{P}, 80 \text{ MHz to } 800 \text{ MHz}$ $d=2,3\sqrt{P}, 800 \text{ MHz to } 2,7 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz	3 V/m	<p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range ^b. Interference may occur in the vicinity of equipment marked with the following symbol:</p>
Proximity fields from RF wireless communications equipment IEC 61000-4-3	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment.</p>			
<p>^b Over the frequency range 150 kHz to 80 kHz, field strengths should be less than 3 V/m.</p>			

INSTALLATION

Risk of incorrect installation.

Only a qualified technician may install the compressor and place it into operation for the first time. Their duty is to train operating personnel on the use and maintenance of the equipment. An entry is made in the equipment installation record to certify installation and operator training. (See warranty card)

10. INSTALLATION CONDITIONS

- The compressor may only be installed and operating in dry, well-ventilated and clean environments under the conditions specified in the Technical Data chapter.

Risk of damage to the device.

The equipment may not be operated outdoors or in otherwise wet or damp environments.

Risk of explosion.

Do not use the equipment in the presence of explosive gases, dust or combustible liquids.

Burn or fire hazard! Caution! Hot surface!

Portions of the compressor, dryer may be hot and reach hazardous temperatures during compressor operation that may harm materials or operating staff.

- The compressor must be installed so that it is accessible at all times for operating and maintenance. Please ensure that the nameplate on the device is readily accessible.
- The compressor must stand on a flat, sufficiently stable base (be aware of the weight of the compressor, see the Technical Data chapter).
- Any other use or use outside this framework is not considered as intended use. The manufacturer is not responsible for any damages arising from such use. The operator/user assumes all risks.

High temperature hazard.

The placement of air flow impediments upstream or downstream of the cooler is prohibited. The temperature of internal and external parts of the cooler may be hot and reach hazardous temperatures.

You may notice a “new product” odour when you first place the product into service (for a short period of time). This odour is temporary and does not impede the normal use of the product. Ensure the space is properly ventilated after installation.

11. COMPRESSOR ASSEMBLY

11.1. Handling and releasing the compressor

- Unpack the compressor from the packaging.
- For compressors with a cabinet, open the front door and remove the connecting strip (26) in the front lower part. Remove the cabinet from the compressor (Fig. 2).
- Position the compressor at the site of installation (Fig. 3).


Use the handles on the compressor to move the compressor as needed. Do not use other parts of the compressor (the air pump, cooler, etc.) for grip

The number of persons required to handle the equipment must be matched to the weight of the device.

Fig. 3: Handling the compressor


- Remove the transport securing elements from the air pumps (Fig. 4).

Prior to installation, ensure that the compressor is free of all transport packaging and stabilizers to avoid any risk of damage to the product.

Remove all devices used to secure the air pumps once the compressor is installed and levelled at the site of final installation.

Fig. 4: Releasing the air pumps

12. PNEUMATIC CONNECTION

12.1. Connecting to the compressed air outlet

- Connect the pressure hose (2) to the compressed air outlet (1) of the compressor.
- Route the hose to the compressed air manifold or directly to the appliance – dental set.
- Secure the pressure hose in the clamp (3) (Fig. 5).
- Route the pressure hose through the opening in the rear wall of the cabinet for cabinet-mounted compressors (Fig. 7).



Fig. 5: Connecting to the compressed air outlet

12.2. Connecting the cabinet pressure gauge to the compressor

- Remove the plug (1) from the threads (2) on the pneumatic block of the compressor.
- Connect the cabinet pressure gauge hose to the threaded fitting.


Fig. 6: Connecting the cabinet pressure gauge to the compressor

12.3. Connecting the condensate vessel

- Connect the condensate outlet of the dryer and the condensate collection vessel with a hose (1). (Fig. 7).
- With the compressor in the cabinet, route the hose through the hole in the back wall of the cabinet and connect it to the collection tank. Mount the condensate collection tank to the side panel or to the back panel of the cabinet (Fig. 7).

Injury hazard.

The condensate drain must not be directly connected to the waste or left without connection to the collection vessel!

Passers-by may be injured!



Fig. 7: Condensate drain


Risk of damage to pneumatic components.

Air hoses must not be broken.

12.4. Floor installation

Connect the compressor using the hoses prepared in advance in the floor, as per the installation plan.

Fig. 8: Floor installation

Description to Fig. 8:

- | | |
|---------------------------|---|
| 1. Outline of the cabinet | 4. Compressed air connection G3/8 |
| 2. Outline of the base | 5. Power supply input
230V/50(60) Hz
3G1.0x4000 |
| 3. Front section - door | |

- All dimensions are in millimetres
- Distance from wall min. 100 mm

13. ELECTRICAL CONNECTION

- The product is delivered with a cord equipped with a plug with earthing pin. Insert the mains plug into a rated mains receptacle.
- Keep the socket easily accessible to ensure that the device can be safely disconnected from the mains.
- Connection to the power distribution cabinet must be max. 16 A.

Risk of electric shock.

It is necessary to follow all local electro technical regulations. The mains voltage and frequency must comply with the data stated on the device label.

Risk of fire and electric shock.

Electrical cord must not be broken.

Risk of fire and electric shock.

Electrical cable must not be in contact with hot compressor components.

- Connect the pin for equipotential bonding 6Ø mm (1) to the electrical circuit using the defined method under valid electrical engineering regulations.
- The socket for equipotential bonding (2) is an accessory and is not provided in the basic package.

Fig. 9: Equipotential bonding

13.1. Connecting a compressor not installed in a cabinet

- Insert the mains plug into a rated mains socket.
- The compressor is ready for operation.


13.2. Connecting a compressor installed in a cabinet

- Route the mains plug through the opening in the rear wall of the cabinet for cabinet-mounted compressors (Fig. 10).
- Connect the cabinet electrically to the compressor by inserting the provided power cord with connector (1) into an outlet (2) (Fig. 11).
- Disconnect the cabinet electrically from the compressor by pulling the connector from the outlet while the latch is released.

Fig. 10: Opening in the rear wall of the cabinet

Fig. 11: Connect the cabinet to the compressor

14. PLACEMENT OF THE COMPRESSOR IN THE CABINET

- Position the compressor within a distance of at least 200 mm from the cabinet to facilitate the movement of the hoses and the power cord within the cabinet (Fig. 12).


Fig. 12: Positioning of the compressor for connections

- Make the pneumatic connections per Chapter 12.1.
- Route the pressure hose, the condensate drain hose, and the power cord through the opening in the rear wall of the cabinet.
- Position the compressor at such a distance from the cabinet so the WINSTA connector can make the electrical connection between the compressor and the cabinet (Fig. 11, Fig. 13).


Fig. 13: Positioning the compressor for electrical connections

- Electrically connect the compressor per Chapter 13.
- Slide the compressor into the cabinet and secure with the connecting strip (26).
- Connect the condensation collection vessel per Chapter **Chyba! Nenašiel sa žiaden zdroj o**

dkazov..

Do not push the compressor all the way to the back of the cabinet as it may permanently damage the dryer.

- Check to ensure the compressor is fully inserted into the cabinet using the control dimension (Fig. 14).



Fig. 14: Positioning of the compressor in the cabinet

14.1. Changing door direction

- Disconnect the earthing conductor from the cabinet. (A)
- Remove the door, the rectifying screws on both sides of the cabinet (3) and the hinge (4) bracket (2). (B)
- Install the bracket for hinge (2) on the left side of the cabinet. (C)
- Rotate the door 180°. (B)
- Insert the spacer (5) between hinge (1) and the bottom of the door. (C)
- Install the door, the 2 rectification screws (3) and connect the earthing conductor. (C)
- Remove the door lock (6) and rotate 180°. (D)
- Remove the latch (7) and rotate 180°. (D)
- Install the lock. (D)

A

15. COMMISSIONING

- Make sure all transport stabilizers were removed.
- Check that all compressed air hose connections are correct (see chap. 12).
- Check correct connection to the mains (see chap. 13).
- Check the circuit breaker position, it must be in the “I” position. If the main switch (4) is in the “0” position, turn it to the “I” position (Fig. 15).
- For cabinet-mounted compressors, turn the main switch (5) on the front side of the cabinet into

the "I" position; a green indicator indicates the equipment is in operation (Fig. 15).

- Check connection of the cabinet connector to the compressor (Fig. 11).
- Check connection of the cabinet manometer hose to the pneumatic block of the compressor (Fig. 6).
- Check the correct connection of the hose from the condensate outlet on the dryer to the collection vessel. (Fig. 7).


The compressor is not equipped with a backup power supply.

16. PNEUMATIC AND ELECTRICAL DIAGRAMS

16.1. Pneumatic diagram

DK50 PLUS/M, DK50 2V/M

Description to pneumatic diagram:

- | | |
|--------------------|----------------------------------|
| 1 Inlet filter | 9 Pressure gauge |
| 2 Air pump | 10 Cooler |
| 3 Fan | 11 Cooler fan |
| 4 Solenoid valve | 12 Condensate separator |
| 5 Safety valve | 13 Dryer |
| 6 Non-return valve | 14 Condensate collection vessel |
| 7 Pressure switch | 15 Regulator with pressure gauge |
| 8 Air tank | |


16.2. Electrical diagrams

DK50 PLUS/M, DK50 PLUS S/M DK50 2V/M, DK50 2V S/M

1/N/PE 230 V, 50/60 Hz

115 V, 60 Hz


ELECTRIC OBJECT OF 1st. CAT.

Compressor cabinet

1/N/PE ~ 230V 50/60Hz

ELECTRIC OBJECT OF 1st. CAT.

Description to electrical diagrams:

M1	Compressor motor	B1	Pressure switch
E1	Compressor fan	X1	Terminal box
E2	Cooler fan	F1	Circuit breaker switch
E10	Cabinet fan	S10	Switch
B2	Temperature switch	P11	Indicator
X10,X11	Connector	M11	Solenoid valve
X2	Connector (Molex)		

Note:

- ** - Connect the jumper only for compressor models not installed in cabinets (see chap. 22.9)
- * - Design only applies to the cabinet version of DK50 PLUS S/M, DK50 2V S/M

OPERATION

ONLY TRAINED PERSONNEL MAY OPERATE THE EQUIPMENT!

Risk of electric shock.

In case of emergency, disconnect the compressor from the mains (pull out the mains plug).

Risk of burn or fire.

During compressor operation, the pump parts may heat to temperatures dangerous for contact with persons or materials.

Warning – compressor is controlled automatically.

Automatic start-up. When pressure in the pressure tank drops to switch-on pressure, the compressor automatically switches on. The compressor automatically switches off, when pressure in the air tank reaches the switch-off pressure.

Hazard of damage to pneumatic components.

Working pressures of the pressure switch set by the manufacturer must not be changed. If the compressor operates at higher working pressure than the switch-on pressure, it means that the air consumption is too high (see Chapter Malfunctions).

Hazard of damage to dryer.

When the dryer is operated at ambient temperature that is higher than maximum operating temperature stated in Technical data, the dryer can be damaged.

Required drying performance can only be achieved when the defined operating conditions are followed.

Drying performance decreases and the achieved dew point drops if the dryer is operated at any pressure below the minimum working pressure.

During longer operation of the compressor, the temperature in the surrounding of the compressor increases above 40 °C and the cooling fan switches on automatically. After cooling the space below 32 °C, the fan switches off.

17. SWITCHING ON THE COMPRESSOR

Start the compressor (without a cabinet) at the pressure switch (1) by turning the switch (2) to position "I." This starts the compressor and fills the tank to the switching off pressure, which then shuts off the compressor.


Start the compressor (with a cabinet) at the switch (5) on the front side of the cabinet, the indicator will turn green. This starts the compressor and fills the tank to the switching off pressure, which then shuts off the compressor.

The compressor operates in automatic mode, switched on and off by the pressure switch, depending on compressed air consumption.

Check the values of switch-on and switch-off pressure using the pressure gauge (3). Tolerance of $\pm 10\%$ is acceptable. Pressure in the air tank must not exceed the allowed operating pressure.

After turning on the compressor, the dryer's program version and then the number of working hours will be displayed for a short time on the display of the dryer's control electronics.

Fig. 15: Switching the compressor

Risk of damage to pneumatic parts.

The pressure switch (1) was set up by the manufacturer and any other changes to its settings may be performed only by a qualified technician trained by the manufacturer.

Compressor – when first started and placed into service, the compressor fills the air tank until the switching pressure and the pressure switch turn off the compressor. The compressor operates in automatic mode, switched on and off by the pressure switch, depending on compressed air usage.

When the compressor is running, the compressed air passes through a cooler and a dryer that removes moisture from the compressed air.

18. SWITCHING OFF THE COMPRESSOR

- Switching off the compressor due to service or any other reason shall be performed using the pressure switch (1) by turning the switch (2) to position „0“ (Fig. 15) and **pulling out the mains plug from the socket**. This disconnects the compressor from the mains supply.
- Vent the pressure in the air tank to zero by opening the drain valve (7). (Fig. 1)

ADS DRYER

19. PRINCIPLE OF OPERATION

The dryer is controlled by a signal from the compressor pressure switch.

The AD dryer for compressors up to 140 l/min (5 bar) is a single-chamber adsorption dryer based on the PSA method ("Pressure Swing Adsorption"). The principle of the PSA method consists in switching the modes of the dryer chamber. The given dryer can be operated for an extended time in the S3 operating mode and for a short time in the S1 mode.

Moist air entering the chamber in the S1 operating mode is dried by means of adsorption in the dryer chamber. After reaching a certain time, the regeneration valve in the lower part of the dryer opens and the chamber briefly switches to the regeneration phase. During this phase, the dried air from the air tank passes through the regeneration nozzle and the adsorbent in the chamber is regenerated, while the air from the compressor aggregate is discharged through the regeneration valve.

Risk of damage to the dryer.

The adsorption dryer is not designed to work in the S1 operating mode!

19.1. Dryer regeneration


Regeneration of the dryer takes place during the break in the compressor operation using the regeneration nozzle.

The nozzle is physically designed for the selected level of regeneration at the maximum considered mode of operation, e.g. S3 50%, S3 70% and so on.



20. BASIC SCREENS

All segments of the display light up for inspection once the power is connected.

The display briefly shows the inscription rSA05 and then the program version.

The output for the solenoid valve is switched on (approx. 5 s) to ensure a smoother start of the compressor motor. Subsequently, the message "M_H" will appear on the display for about 1 s, followed by a number representing the total number of motor hours of the dryer.

Displaying the dryer's total motor hours

After pressing the button, the message "S_H" is displayed for approx. 1 s, followed by a number representing the number of motor hours of the dryer since the last service. The last decimal point is illuminated to distinguish it from the total motor hours.


Displaying the dryer's motor hours since the last service

Once pressed again, the display shows the dryer's total motor hours. The status of the unit is indicated by two LEDs.

The green "STATUS" LED is lit when the compressor motor is running and flashes with a frequency of 1 Hz when the regeneration valve is activated.

The orange "SERVICE" LED flashes with a frequency of 1 Hz when there are less than 100 hours left until service is required.

At the same time, the message "S_A" appears on the display once every 10 s, followed by the number of hours until the service is performed.

Display of remaining motor hours until servicing

The service interval of the dryer is set to 4,000 motor hours, this means that this warning will appear after 3,900 motor hours have passed since the last service and the service hours have been reset. After exceeding the interval of 4,000 motor hours since the last service, the orange "SERVICE" LED remains on, and the display shows "S_O" once every 10 s, followed by the number of hours by which the service interval has been exceeded.

Display of hours past the defined service interval

The unit measures the running time of the compressor motor using a signal from the pressure sensor, from which it then calculates the necessary regeneration time to activate the regeneration valve. If the compressor runs for longer than 3,600 s, the unit starts forced regeneration. This consists of switching on the regeneration valve for 7 s during every 70 s when the compressor is running. The unit also continuously calculates the regeneration deficit, which is stored in the processor's EEPROM memory every 300 s. This deficit is applied as the necessary regeneration time after the compressor is turned off and then turned on.

21. SERVICE SCREENS

For a better overview of the measured parameters and possible diagnostics, the RSA05 unit can be switched to the service parameters display mode. Do this by pressing the button and holding it down for approx. 3 to 5 s.

The service display is indicated by the "t_r" indicator for approx. 1 s and then by the display of the regeneration time.

Display of regeneration time "t_r"

In regeneration mode, the remaining regeneration time in seconds is displayed here, which gradually decreases. With the compressor in operating mode, the remaining time from the last regeneration is displayed.

Briefly pressing on the button scrolls through the individual parameters in the following order.

Display of compressor run time "t_b"

With the compressor in operating mode, the compressor run time in seconds is displayed here, which increases gradually. In the regeneration mode, the last compressor run time is displayed.

Display of regeneration deficit "t_d"

This shows the regeneration deficit in seconds. In regeneration mode, this figure is the same as the remaining regeneration time. With the compressor in operating mode, the regeneration deficit gradually increases as the compressor run time increases. This data is stored in the EEPROM memory every 300 s and is applied once the dryer is activated.

Display of forced regeneration time "t_n"

Here, the time of forced regeneration is displayed in seconds, which starts to increase after the compressor is switched to continuous run mode ($t_b \geq 3600$ s).

Display of total motor hours "M_H"

Here, as in the operating mode, the total running hours of the compressor are displayed.

Display of service hours "S_H"

This display shows the motor hours of the compressor since the last service. After the service has been performed, the service hours can be reset to zero in the service screen. Press and hold the button for approx. 3 to 5 s to reset.

Display of motor minutes "M_M"

The compressor motor minutes are shown here and are stored in the EEPROM memory every 30 minutes.

The service screen automatically closes after approx. 30 minutes have passed since the last button was pressed. The unit will go to the standard operating screen showing the compressor motor hours.

PRODUCT MAINTENANCE

22. PRODUCT MAINTENANCE

The operator shall ensure completion of repeated testing of the device at least once every 24 months (EN 62353) or at intervals defined by applicable national legal regulations. A record of these test results shall be completed (e.g. per EN 62353, Annex G) together with the measurement methods.

The equipment has been designed and manufactured to keep maintenance to a minimum. The following work must be performed to preserve the proper and reliable operation of the compressor.

Unauthorised interference hazard.

Repair work outside the framework of standard maintenance (see Chapter 22.1) may only be performed by a qualified technician (an organisation authorized by the manufacturer) or the manufacturer's customer service.

Standard maintenance work (see Chapter 22.1) may only be performed by the operator's trained personnel.

Only use manufacturer-approved replacement parts and accessories.

Danger of injury or equipment damage.

Prior to commencing compressor maintenance, it is necessary to:

- check if it is possible to disconnect the compressor from the appliance in order to avoid any risk of injury to the person using the appliance or other material damage;
- turn off the compressor;
- disconnect it from the mains (pulling the cord out of the mains socket);
- vent the compressed air from the air tank.

Venting compressed air poses an injury hazard.

Wear eye protection, i.e. goggles, when venting compressed air from the compressed air circuit (air tank).

Burn hazard.

When the compressor is running or shortly thereafter, certain portions of the air pump may be hot - do not touch these components.

Injury hazard.

Allow the equipment to cool down before maintenance, service or connecting/disconnecting the compressed air supply!

The removed grounding conductor during service must be connected back to the original position after completing the service.

- For any maintenance or repair work, the compressor may be removed from the cabinet by pulling it from the cabinet (on rolling casters as much as the cord between the compressor and the cabinet allows). (see chap. 22) The required maintenance or repair work may then be performed.
- If the cord between the compressor and the cabinet becomes disconnected (by pulling the mains plug, see chap. 13.2), the compressor shall not be functional. In order to bring the compressor into function again, please follow, see chap. 13.2 or 22.9 or 22.10.

22.1. Maintenance intervals

Time interval ^{b)}	50 Hz	60 Hz	Once a day	Once a week	Once a year	Once every 2 years	Time intervals (h)						Chap.	Set of replacement parts	Performed by
							1600 h	2000 h	3200 h	4800 h	6400 h	8000 h			
Check of product operation	x														operator
Drain the condensate from the air tank - At high RH) - Pour the condensate out of the vessel	x														operator
Drain the condensate from the air tank - At normal RH		x													operator
Product function check			x												operator
Leak check connections and inspect device					x										qualified technician
Inspection of electrical connections					x										qualified technician
Cooler and fan check					x										qualified technician
Safety valve check					x										qualified technician
Replacing the inlet filter on the DK50 PLUS air pump ^{a)}					x				x					025200126-000	qualified technician
Replacing the inlet filter on the DK50 2V air pump ^{a)}									x					025200139-000 025200150-000	qualified technician
Replacement of cassette with adsorbent media ADS 70 P dryer (DK50 PLUS/M)										x				603032090-000	qualified technician
Replacement of cassette with adsorbent media ADS 140 dryer (DK50 2V/M)											x			603022655-000	qualified technician
Conduct a " repeated test" per EN 62353						x									qualified technician

^{a)} data is in hours, if not available, data is in years

^{b)} time interval reduces by 20 % for compressor variants 60 Hz (2000h./1600h., 4000h./3200h., 6000h./4800h., 8000h./6400h., 10000h./8000h., 12000h./9600h.)

22.2. Check of product operation

- Check air pump condition – the air pump should be operating normally without excessive vibration or noise. Troubleshoot any problem or call in service personnel if trouble is detected.
- Visually inspect fan operation – the fans must be operating when the aggregates are running. Troubleshoot any problem or call in service personnel if trouble is detected.
- Check to ensure the power cord, the connecting compressed air hoses are undamaged. Replace damaged components or call in service personnel.
- Check the ambient temperature – the ambient temperature must be below the temperature limit (40°C). Cool the space if the temperature is high.
- Compressor with air dryer - open the plug on the condensate collecting vessel and drain the condensate.
- Check the operating condition of the equipment (see chap. 22.4).

22.3. Check the compressed air connections for leaks and inspect the equipment

Leak testing

- Check the compressor's compressed air lines for leaks during operation – pressure supplied by the compressor.
- Use a leak analyser or soapy water to check all joints and connections for leaks. Tighten or reseal the connection where leaks are found.

Inspecting the equipment

- Check the condition of the compressor aggregate for normal operation and noise levels.
- Fan operation check – the fans must be running during the defined compressor work cycles.
- Check function of the temperature switch (B2) – heat the temperature switch to temperature above 40°C (e.g. using a heat gun – Be careful, do not operate with high temperature on plastic components in the surrounding, it can cause deformation of plastic materials). Fan E1 - the compressor must be under voltage.
- Check filters condition – filters must be undamaged and sufficiently clean.
- Check condition of the pump itself, check if there is no dirt inside the crankcase or clearance in the crankshaft.
- Check the operation of the automatic condensate drain.
- Replace any defective parts as needed.
- Check the operating condition of the equipment (see Chapter 22.4).

22.4. Inspection of electrical connections

Risk of electric shock.

Inspect the product's electrical connections when the mains are disconnected.

- Check mechanical function of the main switch.
- Check if the power cable and conductors are not damaged.
- Visually check if cables are connected to the terminal box.
- Check all screw connections of the green-yellow PE grounding conductor.

22.5. Condensate drain

A wet floor resulting from overflow from the vessel poses a slip hazard.

Regularly empty the condensate collection vessel.

Injury hazard.

Do not connect the condensate drain directly to a waste drain!

Passers-by may be injured.

Switch off the compressor every time before emptying the vessel!

Condensate is automatically separated into the condensate collection tank.

- Monitor the level in the vessel using the markings (depending on the volume of the vessel), and empty at least once a day.

Fig. 16: Check of the condensate collection vessel

Injury hazard.

Carefully and slowly open the drain valve. Quickly opening the drain valve is accompanied by a loud sound and an uncontrollable ejection of the accumulated condensate.

If necessary, it is possible to drain the air tank using the drain valve (1) from the lower side of the air tank (Fig. 17).

Switch off the compressor at the mains and reduce the air pressure in the device to a maximum of 1 bar, e.g. by venting air through the connected equipment.

Use a hose to drain with the free end placed in a suitable container and then open the valve to drain condensate from the tank.


Fig. 17: Condensate drain

Before the following checks it is required:

- **Compressor with cabinet** – release the door lock and open the cabinet door.

22.6. Check of safety valve

- Turn the screw (2) on the safety valve (1) several times to the left until the safety valve releases the air.
- Let the safety valve vent for a few seconds.
- Turn the screw (2) fully to the right, the valve must be closed now.


Fig. 18: Check of safety valve

Damage to the safety valve could cause pressure to rise to hazardous levels.

Never use the safety valve to release the air pressure in the air tank. This could damage the safety valve. The valve is set to the maximum permitted pressure by the manufacturer.

Never adjust a safety valve.

Venting compressed air poses an injury hazard.

Wear safety glasses when inspecting a safety valve.

22.7. Inlet filter replacement

The inlet filter (1) is located inside the compressor crankcase cover.

Inlet filter replacement:

- Pull out the rubber plug (2) by hand.
- Remove the dirty inlet filter (1).
- Insert a new filter and mount the rubber plug.


Fig. 19: Inlet filter replacement

Inlet pre-filter replacement

- Pull out the pre-filter by hand (3).
- Replace the element and re-insert.


Fig. 20: Inlet pre-filter replacement

22.8. Inspecting the cooled and fan

The equipment, in particular the compressor fan, cooler fan, and the cooler, must be kept clean to ensure efficient drying (Fig. 1) - remove dust from the surface of the cooling fins and fans by vacuuming or blowing down with compressed air.

22.9. Procedure for connecting a compressor disconnected from the cabinet

Prior to any maintenance or repair work, switch off the compressor and disconnect it from the mains (pull out the mains plug).

The compressor (without cabinet) requires that the jumper is always mounted in the terminal strip for correct operation (Fig. 21, pos. B). It then replaces the switch on the cabinet.

If there is no jumper in the terminal strip of the compressor (disconnected from the cabinet), the compressor will not work!

If the cord between the compressor and the cabinet becomes disconnected (by pulling the mains plug) and the compressor is removed from the cabinet, the compressor will not be functional. This is why it is necessary to first create a connection to the terminal strip with a jumper (replacing the circuit breaker function) (Fig. 21).



PROCESS:

Jumper installation (A-B process)

- Disconnect the product from the mains by pulling the plug out of the socket.
- Remove the cover on the electrical panel (on the compressor).
- The jumper is not in the terminal strip – A.
- Insert the jumper into the terminal strip – B.
- Reinstall the cover on the electrical panel.
- Connect the device to the electrical mains.
- Activate the compressor by turning on the switch on the pressure switch.

Compressor - when operating outside the cabinet, the electrical connection with an installed jumper **must** be created (Fig. 21, pos. B)

Fig. 21

230V

22.10. Procedure for connecting a compressor to a new cabinet

Prior to any maintenance or repair work, switch off the compressor and disconnect it from the mains (pull out the mains plug).

The compressor in a cabinet requires that the jumper is not mounted in the terminal strip for correct operation (Fig. 21, pos. A). The switch on the cabinet switches the entire product on and off, including the compressor.

If there is a jumper in the terminal strip of the compressor connected to a cabinet, the switch on the compressor cabinet will not work!

To connect a new cabinet to a compressor that previously worked without a cabinet, it is necessary to cancel the connection in the compressor terminal strip by performing the procedure as follows. (Fig. 21). See also chap. 13.1 and 13.2.

PROCESS:

Jumper removal (B-A process):

- Disconnect the product from the mains by pulling the plug out of the socket.
- Remove the cover on the electrical panel (on the compressor).
- The jumper is in the terminal strip – B.
- Remove the jumper from the terminal strip – A.
- Reinstall the cover on the electrical panel.
- Insert the compressor into the cabinet.
- Connect the cabinet to the compressor by inserting the power cord with connector (1) into a power socket (2) (Fig. 11).
- Connect the device to the electrical mains.
- Activate the compressor by turning on the switch on (2) the pressure switch (1) and the switch (5) on the cabinet (Fig. 15).

Compressor with cabinet - when operating in the cabinet, there **must not** be an electrical connection (Fig. 21, pos. A).

Note: **In case of using a compressor without a cabinet or in case of maintenance, a jumper which is not installed into the terminal strip should be kept** in a designated place next to the electrical panel.

22.11. Cleaning and disinfection of the exterior surfaces of the product

Clean and disinfect the exterior surfaces with neutral cleaning products.

Use of aggressive detergents and disinfectants containing alcohol and chlorides can lead to surface damage and discolouration.

23. LONG-TERM SHUTDOWN

If the compressor is not going to be used for long period, it is recommended to drain all condensate from the air tank and switch on the compressor for about 10 minutes, keeping the drain valve open (1). Thereafter switch off the compressor using the switch (2) on the pressure switch (1) (Fig. 15), close the drain valve and disconnect the device from the mains.

24. DISPOSAL OF DEVICE

- Disconnect the equipment from the mains.
- Release the air pressure in the pressure tank by opening the drain valve.
- Dispose of the equipment following all applicable regulations.
- Entrust a specialised company to sort and dispose of waste.
- Worn out components have no negative environmental impact

TROUBLESHOOTING**Risk of electric shock.****Before interfering with the equipment, first disconnect it from the mains (remove the power socket).****Working with pressurised pneumatic components poses a risk of injury.****Before interfering with the equipment, vent the air tank and the compressed air system to zero pressure.****Troubleshooting may only be performed by a qualified service technician.**

Malfunction	Possible cause	Solution
Compressor does not switch on	No voltage in the pressure switch	Check voltage in socket
		Check circuit breaker switch – switch to position switched-on „I“
		Loose conductor from terminal - repair
		Check electrical cord – replace defective cord
	Motor winding failure, damaged thermal protection	Replace motor or windings
Capacitor failure	Replace capacitor	
Seized piston or other rotating part	Replace damaged components	
Pressure switch does not switch	Check function of pressure switch	
The cabinet switch does not switch off the compressor	There is a compressor in the cabinet that previously worked separately (there is a jumper in the compressor terminal strip)	Remove the jumper from the terminal strip, see chap. 22.10.
Compressor often switches on	Air leakage in pneumatic system	Check pneumatic system – seal loose connections
	Non-return valve (SV) leakage	Clean non-return valve, replace seals, replace non-return valve
	Large volume of condensed liquid in pressure tank	Drain condensed liquid
	Low compressor capacity	Check time of filling the air tank
Low pressure in the air tank (compressor running constantly)	High air consumption of supplied equipment	Decrease air consumption Use compressor with higher capacity
	Leakage in pneumatic system	Check pneumatic system – seal loose connections
	Low capacity of the pump	Clean / replace the pump
	Pump failure	Clean / replace the pump
	Dryer failure	Replace dryer

Prolonged operation of the compressor	Air leakage in pneumatic system	Check pneumatic system – seal loose connection
	Worn piston ring	Replace worn piston ring
	Dirty inlet filter	Replace dirty filter with new filter
	Incorrect function of solenoid valve	Repair or replace fan or coil
Compressor is noisy (knocking, metal noises)	Damaged piston bearing, piston rods, motor bearing	Replace damaged bearing
	Loose (cracked) dampening element (spring)	Replace damaged spring
Dryer does not dry (condensed water in the air) *	Defective cooler fan	Replace fan Check power source
	Damaged dryer	Replace dryer
Degraded drying performance – high-pressure dew point (condensed water in the air)	Low operating pressure	Reduce the demand for air, check the output from the compressed air source, fix any leaks in the distribution system
	Regeneration solenoid valve not working	Check coil operation, replace if damaged
		Inspect the condition of the valve – clean the valve or replace if problems persist
	Air regeneration nozzle plugged	Clean or replace the nozzle (see product maintenance)
	Cooling fan not working	Check the power source to the fan Replace damaged fan
	Dirty cooler	Inspect the cooler and clean as necessary
	Clogged silencer in the condensate collection tank	Check the condition of the silencer. In the case of high flow resistance, replace the silencer
	Defective dryer electronics	Check the display electronics, and if the display shows the compressor run time, check the connection of the control electronics. In case of malfunction, replace the electronics
The compressor operated in a higher mode (S3 > 80%)	Keep the compressor running in a mode that is less than S3-50%. The dryer will automatically regenerate. Do not switch off the compressor!	
Dryer emitting high levels of noise	Damaged fan	Replace damaged fan
	Damaged silencer	Replace the silencer
	High pressure in equipment due to excessively clogged filters	Change the cartridges in the dryer

)* Clean inside surfaces of the air tank thoroughly and remove all condensed liquid after dryer failure.

Once a dryer fault is cleared and after its reassembly, the dryer should quickly be regenerated, best when using continuous compressor operation at a pressure of around 1 bar for a period of at least 1 hour; then check the dryness of the compressed air.

Check the moisture content of the air exiting the air tank (see the Technical data chapter) to prevent damage to connected downstream equipment.

25. REPAIR SERVICE

Warranty and post-warranty repairs must be done by the manufacturer, its authorized representative, or service personnel approved by the supplier.

Attention.

The manufacturer reserves the right to make changes to the equipment without notice. Any changes made will not affect the functional properties of the equipment.

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- Držiak pántu (2) namontovať na ľavú stranu skrinky. ©