DK50 PLUS MOBILE DK50 2V MOBILE

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OBSAH 71

CONTENTS

GENERAL INFORMATION.....	6
1. CONFORMITY WITH THE REQUIREMENTS OF THE EUROPEAN UNION.....	6
2. SYMBOLS.....	6
3. DEVICE USE	7
4. GENERAL SAFETY INSTRUCTION	8
5. STORAGE AND TRANSPORT	9
PRODUCT DESCRIPTION	10
6. VARIANTS	10
7. PRODUCT FUNCTION	10
TECHNICAL DATA.....	14
INSTALLATION.....	18
8. INSTALLATION CONDITIONS	18
9. COMPRESSOR ASSEMBLY	18
10. PNEUMATIC CONNECTION	19
11. ELECTRICAL CONNECTION	19
12. COMMISSIONING	20
13. PNEUMATIC AND ELECTRICAL DIAGRAMS	21
OPERATION.....	27
14. SWITCHING ON THE COMPRESSOR.....	28
15. SWITCHING OFF THE COMPRESSOR.....	28
PRODUCT MAINTENANCE	29
16. PRODUCT MAINTENANCE.....	29
17. LONG-TERM SHUTDOWN.....	35
18. DISPOSAL OF DEVICE	35
TROUBLESHOOTING	36
19. REPAIR SERVICE	37
ANNEX	104
20. INSTALLATION RECORD	104

GENERAL INFORMATION

Read the User manual carefully and keep it before use of the product. The User manual provides information on correct use – installation, operation and maintenance of the product.

The User manual corresponds with the design of the product and condition according to the applicable safety and technical standards at the time of printing. The manufacturer reserves all rights concerning the protection for the stated connections, procedures and names.

Slovak version is original of the User manual. The User manual has been translated in accordance with the best available knowledge. The Slovak version is to be used in the event of any uncertainties.

This user manual is the original instructions. Translation is performed in accordance with the best available knowledge.

1. CONFORMITY WITH THE REQUIREMENTS OF THE EUROPEAN UNION

This product conforms to the requirements of the European Union 2006/42/EC, 2014/29/EU, 2014/35/EU, 2014/30/EU, 2011/65/EU and is safe if used in compliance with the intended use and if all safety instructions are followed.

User manual is in compliance with requirements of Directive 2006/42/EC.

2. SYMBOLS

The following symbols and marks are used in the User manual, on the device and its packaging:

General warning

Warning – risk of electric shock

Warning - compressor is controlled automatically

Warning – hot surface

General caution

Refer to instruction manual

CE – marking

Serial number

Article number

Protecting earthing

Terminal for ground connection

	Compressed air outlet
	Package handling label – fragile
	Package handling label – this side up
	Package handling label – keep dry
	Package handling label – temperature limits
	Package handling label – limited stacking
	Package label – recyclable material
	Manufacturer

3. DEVICE USE

3.1. Intended use

The compressor is used as source of clean oil-free compressed air intended to be used in industry and laboratories, where parameters and properties of the compressed air are suitable.

The compressor is exclusively intended to compress air without content of explosive or chemically unstable substances.

The compressor is intended for operation in clean and dry rooms

3.2. Incorrect use

Contamination risk.

Air from the compressor is without additional treatment not suitable for breathing and direct contact with food

Explosion risk.

The product is not intended for operation in rooms with explosion risk.

The compressor must not be used to compress aggressive gases.

The compressor must not be operated in premises with occurrence of flammable vapors.

The compressor must not be operated in other conditions as mentioned in Technical data.

Any other use of the product beyond the intended use is considered as incorrect use. The manufacturer is not responsible for any damages or injuries as a result of incorrect use or disobedience to instructions stated in this User manual. All risks shall be solely borne by the user/operator.

4. GENERAL SAFETY INSTRUCTION

The product is designed and manufactured so that any risks connected with its use are minimized and the product is safe for the user and surrounding when used according to the intended use and the instructions stated below are followed.

4.1. Required qualification of the personnel

- Each user must be trained by the manufacturer or an organization authorized by the manufacturer or instructed on the device operation by other trained user.
- Installation, new settings, changes, extensions and repairs of the product may be performed by the manufacturer or an organization authorized by the manufacturer (hereinafter qualified technician).
- Otherwise the manufacturer is not responsible for safety, reliability and correct functioning of the product.

4.2. General instruction

- When operating the compressor, all acts and local regulations valid in the place of use must be observed. The operator and user are responsible for following the applicable regulations.
- Before every use, the user must check, if the device is functioning correctly and safely. Before building the compressor in other devices, the supplier must assess, if the supplied air and construction of the device comply with the requirements of the specified intended use. Taking this into account, follow the product technical data. Assessment of conformity shall be performed by the manufacturer – supplier of the final product.

4.3. Protection from dangerous voltage and pressure

- The device can be connected only to a properly mounted socket with protective earthing.
- Before connecting the product, check if mains voltage and mains frequency stated on the product are in compliance with the values of the mains.
- Before putting the product into operation, check eventual damages to connected pneumatic hoses and electrical cables. Replace damaged pneumatic hoses and electrical cables immediately.
- Immediately disconnect the product from the mains (pull out the mains plug) in hazardous situations or technical disorders.
- Safety valve settings must not be changed and used to depressurize the air tank.

4.4. Original spare parts and accessories

- Safety of operating staff and failure-free operation of the product are guaranteed only when original spare parts are used. Only accessories and spare parts stated in the technical documentation or expressly approved by the manufacturer may be used.
- The guarantee does not apply to damages arising from use of accessories and spare parts other than prescribed or recommended by the manufacturer and the manufacture is not responsible for them.

5. STORAGE AND TRANSPORT

The manufacturer ships the compressor in a transport packaging. This protects the device from damage during transport.

Risk of damage to pneumatic parts.

The compressor may be transported only depressurized. Vent air pressure from the pressure tank and pressure hoses and drain condensate from the air tank before transporting the compressor.

Original packaging must be kept for eventual return of the device. If possible, always use the original compressor packaging for optimal protection of the product. If it is necessary to return the product within the guarantee period, the manufacturer does not guarantee for damages caused by incorrect packaging of the product.

Transport the compressor in an upright position, always secured by transport fixing elements.

Protect the compressor from high humidity, dirt and extreme temperatures during transport and storage. Do not store together with volatile chemical substances.

If it is not possible to keep the original packaging, dispose of the packaging in an environmentally friendly way. Transport cardboard can be recycled with old paper.

It is not permitted to store and transport the device outside the defined conditions, see below.

5.1. Ambient conditions

Products can be stored in rooms and means of transport that are free from any traces of volatile chemical substances under the following climatic conditions:

Temperature	from -25°C o $+55^{\circ}\text{C}$, in 24h to $+70^{\circ}\text{C}$
Relative humidity	max. 90% (without condensation/ non-condensing)

PRODUCT DESCRIPTION

6. VARIANTS

The compressor is manufactured according to its intended application in the following variants:

DK50 PLUS MOBILE	a compressor in a cabinet with effective soundproofing
DK50 PLUS/M MOBILE	a compressor in a cabinet with an air dryer
DK50 2V MOBILE	a compressor in a cabinet with effective soundproofing
DK50 2V/M MOBILE	a compressor in a cabinet with an air dryer

DK50 PLUS MOBILE
DK50 PLUS/M MOBILE

DK50 2V MOBILE
DK50 2V/M MOBILE

7. PRODUCT FUNCTION

7.1. Compressor without dryer

Fig. 1 The air pump (1) draws in air through an inlet filter (8) and compresses it through a check valve (3) into an air tank (2) from which the appliance draws compressed air. If the pressure in the air tank drops to the switching pressure, the pressure switch (4) turns on the compressor and the compressor supplies compressed air to the air tank until it reaches the switch off pressure, at which time the compressor switches off. The pressure hose is vented through the relief solenoid valve (12) once the compressor unit is shut off. The safety valve (5) prevents the pressure in the air tank from rising above the maximum permitted value. Condensate from the air tank is automatically drained at regular intervals into the vessel via the condensate drain solenoid valve (13). Compressed, oil-free filtered air is stored in the air tank ready for use

7.2. Compressor with membrane dryer

Fig. 2 The air pump (1) draws in air through an inlet filter (8) and compresses it, feeding it to the cooler (10) through the filter (11) and the check valve (3) into the dryer (9), and feeds the clean and dry air into the air tank (2). A portion of the air bypasses the dryer with the captured moisture, which is manifested as a gentle stream of air along the body of the dryer (9). Condensate from the filter is automatically drained at regular intervals into the vessel via the condensate drain solenoid valve (13). The dryer continuously dries the compressed air. Condensate is drained from the air tank via the drain valve (7), depending on the drying controls. Compressed, oil-free filtered air is stored in the air tank ready for use.

For the membrane dryer to work as intended, the dryer is connected directly to the air tank without

the use of a check valve. As the membrane dryer contains selective membrane fibres, there is a small amount of air leakage through them (about 0.5 bar/30 min)

This is why there is a gradual, slow drop in pressure even when there is no demand for compressed air (from the air tank), however, this is not a malfunction.

The pressure vessel does not need to be drained.

7.3. Compressor cabinet

Fig. 3 The cabinet itself is equipped with casters (2 braked and 2 free-rolling) so that the internal shape of the cabinet defines the position of the compressor itself. The rolling casters and the rigid structure of the cabinet itself ensure easy and simple manipulation of the cabinet. The cabinet is composed of a base (on the casters) and a cover connected by 4 columns and enclosed with 4 panels (doors). Every panel (door) may be removed by unscrewing the 2 quick closures.

The top of the front cover of the cabinet is equipped with a pressure gauge showing output pressure, a pressure regulator, main switch, ON/OFF indicator (green indicator) and high temperature indicator for the cabinet (orange indicator). The power cord is routed outside the bottom of the cabinet in the rear (19). A $\varnothing 6/4$ compressed air hose removes condensate from the dryer with an adapter at the condensate drain from the cabinet (24) and the compressed air outlet from the compressor cabinet (23) is connected by a $\varnothing 10/8$ hose (34) to the pressure regulator (33) on the cabinet. The compressed air output with an internal G1/4 thread (or a quick connector for a $\varnothing 12/9$ hose) is located on the surface of the cabinet (20).

The compressor runs in automatic mode when the power switch (28) on the cabinet is switched on.

The following are located on the rear of the cabinet: cooling fans (25) to exhaust hot air from the cabinet through the hot air exhaust ventilation openings in the cabinet (21). The cooling fans are activated at 40 °C by a switching temperature switch and remain on until the temperature in the cabinet falls below 32 °C. At elevated temperatures (above 80 °C) inside the cabinet, a COOLING FAILURE is indicated by the temperature switch via the orange indicator on the panel (30).


Risk of compressor overheating.

Make sure that there are no obstacles at the cooling air inlet into the cabinet (around the bottom part of the cabinet) and at the hot air outlet on the top back side of the cabinet.

Fig. 1: Compressor without dryer

Description for Figures 1-2

1. Air pump
2. Air tank
3. Non-return valve
4. Pressure switch
5. Safety valve
6. -
7. Drain valve
8. Inlet filter
9. Dryer
10. Dryer cooler
11. Filter
12. Solenoid valve
13. Condensate drain solenoid valve
14. Fan
15. Condensate collection vessel
16. Magnetic holder
17. Compressed air outlet
18. Circuit breaker switch

Fig. 2: Compressor with MD dryer



Fig. 3: Cabinet


Cabinet (without panels) – rear view

- 19. Socket
- 20. Compressed air outlet:
 - a. internal thread G 1/4"
 - b. quick coupling for hose Ø12/9
- 21. Heated air output holes
- 22. Tunnel of ventilator for dryer
- 23. Output of compressed air from compressor via hose with Ø 10/8
- 24. Output from a filters for condensate outlet via hose with Ø 6/4
- 25. Ventilators
- 26. Magnetic vessel holder
- 27. Condensate collection vessel

Cabinet (without panels) with compressor DK50 2V/M

- 28. Mains switch
- 29. Green indicator of mains switch turning on
- 30. Orange indicator of operation mode at raised temperature
- 31. Hour counter
- 32. Pressure gauge
- 33. Pressure regulator
- 34. Place for input of hose with Ø 10/8 into pressure regulator from output of compressed air from compressor

20. Compressed air outlet

a. internal thread G 1/4"

b. quick coupling for hose Ø12/9

TECHNICAL DATA

Compressors are designed to operate in dry, ventilated and indoor dust-free rooms with the following climatic conditions:



Temperature	from +5°C to +40°C
Relative humidity	max. 70%

Working pressure 6 – 8 bar		DK50 PLUS	DK50 PLUS/M
Nominal voltage	V, Hz	230, 50/60	230, 50/60
Frequency ^{a)}			
Performance at 6 bar (FAD)	l/min	70/80	60/70
Working pressure ^{b)}	bar	6.0 – 8.0	6.0 – 8.0
Max. current	A	4.1/5.1	4.2/5.2
Motor power	kW	0.55	0.55
Air tank volume	l	25	25
Air quality - filtration	µm	-	-
Maximum operating pressure of safety valve	bar	11.0	11.0
Noise level at 5 bar (L _{pA})	dB	≤57/≤59	≤57/≤59
Operating mode	%	S1-100	S1-100
PDP drying performance at 7 bar	°C	-	≤ +3°C
Time to fill air tank from 0 to 7 bar	s	149/126	198/159
Dimensions (net) w x d x h	mm	550x720x995	550x720x995
Net weight	kg	100 ^{c)}	105 ^{c)}

a) State the compressor variant in the order

b) For other range of pressure consult with the supplier

c) Weight value is only informative data and applies only to a product without any accessories

Dependence of compressor capacity on working pressure

Working pressure 6 – 8 bar		DK50 2V		DK50 2V/M	
Nominal voltage Frequency ^{a)}	V, Hz	230, 50	115, 60	230, 50/60	115, 60
Performance at 6 bar (FAD)	l/min	135		104	
Working pressure ^{b)}	bar	6.0 – 8.0		6.0 – 8.0	
Max. current	A	8.8	20.1	8.9/9.6	20.2
Motor power	kW	1.2	1.1	1.2	1.1
Air tank volume	l	25		25	
Air quality - filtration	µm	-		0.1	
Maximum operating pressure of safety valve	bar	11.0		11.0	
Noise level at 5 bar (L _{pA})	dB	≤62	≤64	≤62	≤64
Operating mode	%	S1-100 S3-60 ^{d)}		S1-100 S3-60 ^{d)}	
PDP drying performance at 7 bar	°C	-		≤ +3°C	
Time to fill air tank from 0 to 7 bar	s	61		78	
Dimensions (net) w x d x h	mm	550x720x995		550x720x995	
Net weight	kg	104 ^{c)}		110 ^{c)}	



^{a)} State the compressor variant in the order

^{b)} For other range of pressure consult with the supplier

^{c)} Weight value is only informative data and applies only to a product without any accessories

^{d)} Applies to voltage version of 115V, 60Hz

Dependence of compressor capacity on working pressure


Working pressure 8 – 10 bar		DK50 2V
Nominal voltage Frequency ^{a)}	V, Hz	230, 50
Performance at 8 bar (FAD)	l/min	110
Working pressure ^{b)}	bar	8.0 – 10.0
Max. current	A	86
Motor power	kW	1.2
Air tank volume	l	25
Air quality - filtration	µm	-
Maximum operating pressure of safety valve	bar	11.0
Noise level at 5 bar (L _{pA})	dB	≤62
Operating mode	%	S1-100
PDP – drying performance at 7 bar	°C	-
Tie to fill air tank from 0 to 7 bar	s	96
Dimensions (net) w x d x h	mm	550x720x995
Net weight	kg	107 ^{c)}

^{a)} State the compressor variant in the order

^{b)} For other range of pressure consult with the supplier

^{c)} Weight value is only informative data and applies only to a product without any accessories

Dependence of compressor capacity on working pressure

FAD correction of capacity for altitude

Capacity given in the form of FAD („Free Air Delivery“) applies to the following conditions:

Altitude	0 m.n.m.	Temperature	20°C
Atmospheric pressure	101325 Pa	Relative humidity	0%

To calculate FAD compressor capacity in dependence on altitude, it is necessary to apply correction factor according to the following table:

Altitude [m.n.m.]	0 -1500	1501 - 2500	2501 - 3500	3501 - 4500
FAD correction factor	1	0,8	0,71	0,60

INSTALLATION

Risk of incorrect installation.

Only a qualified professional can install the compressor and place it into operation for the first time. His obligation is to train the operating personnel on the use and maintenance of the device. He shall confirm installation and training of operators by an entry into the installation record (see Annex).

8. INSTALLATION CONDITIONS

- The compressor may only be installed and operated in dry, well-ventilated and dust-free rooms with conditions stated in Technical data.

Risk of damage to the device.

Do not operate the device in outdoor environment or in humid or wet rooms.

Explosion risk.

Do not install the device in any rooms, in which explosive gases, dusts or flammable liquids may be present.

- The compressor must be installed in such a way, that it is accessible at all times for operation and maintenance and label on the device is accessible.
- The compressor must stand on a flat, sufficiently stable base (be aware of the compressor weight, see Technical data).

You may notice a “new product” odour when you first place the product into service (for a short period of time). This odour is temporary and does not impede the normal use of the product. Ensure the space is properly ventilated after installation.

9. COMPRESSOR ASSEMBLY


9.1. Manipulation and fixation removal

- Unpack the compressor from the packaging.
- Put the compressor at the place of operation.
- Remove transport fixing elements from the pumps (Fig. 4).

Remove fixing elements of the pumps once the compressor has been assembled and balanced at the final installation place.

Fig. 4: Releasing the pump

10. PNEUMATIC CONNECTION

10.1. Connecting to the compressed air outlet (Fig. 3)

There is a G1/4" outlet with internal thread (20) (or quick connection for Ø12 / 9 hose) for the connection of pressure hose in the rear upper part of the housing.

Risk of damage to pneumatic components.

Air hoses must not be broken.

11. ELECTRICAL CONNECTION

- The product is delivered with a cord equipped with a plug with earthing pin.
- Keep the socket easily accessible to ensure that the device can be safely disconnected from the mains.
- The corresponding circuit must be protected in the power distribution max. 16 A.

Risk of electric shock.

it is necessary to follow all local electrotechnical regulations. The mains voltage and frequency must comply with the data stated on the device label.

Risk of fire and electric shock.

Electrical cable must not be in contact with hot compressor components.

Risk of fire and electric shock.

Electrical cord must not be broken.

12. COMMISSIONING

- Make sure all transport stabilizers were removed.
- Check that all pressurized air line connections are secure.
- Check to ensure power is connected correctly to the compressor.
- Turn the switch (28) on the front of the cabinet, to the "I" - green LED indicates the status of the device in operation. (Fig. 3)


The compressor is not equipped with a backup power supply.


13. PNEUMATIC AND ELECTRICAL DIAGRAMS

13.1. Pneumatic diagram

DK50 PLUS, DK50 2V

DK50 PLUS/M, DK50 2V/M

Description to pneumatic diagrams:

- | | |
|--------------------|---------------------------------------|
| 1 Inlet filter | 10 Air tank |
| 2 Compressor | 11 Drain valve |
| 3 Fan | 12 Cooler |
| 4 Relief valve | 13 Coalescing filter |
| 5 Noise dampener | 14 Membrane dryer |
| 6 Non-return valve | 15 Condensate collection vessel |
| 7 - | 16 Solenoid valve of condensate drain |
| 8 Pressure switch | 17 Outlet valve |
| 9 Safety valve | 18 Cabinet pressure gauge |


13.2. Electrical diagrams

DK50 PLUS MOBILE / AOK

DK50 2V MOBILE / AOK


1/N/PE 230 V, 50/60 Hz

ELEKTRIC OBJECT OF 1st. CAT.

DK50 PLUS MOBILE / MD

DK50 2V MOBILE / MD


DK50 PLUS MOBILE / MD

DK50 2V MOBILE / MD

1/N/PE 230 V, 50/60 Hz


ELEKTRIC OBJECT OF 1st. CAT.


DK50 PLUS MOBILE / AOK

1/N/PE 115 V, 60 Hz


ELEKTRIC OBJECT OF 1st. CAT.


DK50 PLUS MOBILE / MD


DK50 2V MOBILE / AOK

DK50 2V MOBILE / MD

Cabinet MOBILE

Description to electrical diagrams:

M1	Compressor motor	F1	Breaker switch
E1	Compressor fan	F11	Breaker
E2	Dryer fan	Q11	Contactator
M10	Relief valve	P11	Hour counter
B1	Pressure switch	Q1	Switch
B2	Temperature switch	P1,P2	Indicator
M4	Condensate drain valve	E10,E11	Cabinet fan
X1	Terminal strip	X10,X11	Connector
X2,X3	Connector		
C1,C _{B1}	Capacitor		

OPERATION

THE EQUIPMENT MAY ONLY BE OPERATED BY THE TRAINED STAFF!

Risk of electric shock

In case of any danger, disconnect the compressor from the mains (pull out the mains plug).

Risk of burn or fire.

During compressor operation, the pump parts may heat to temperatures dangerous for contact with persons or materials.

Warning – compressor is controlled automatically.

Automatic start-up. When pressure in the pressure tank drops to switch-on pressure, the compressor automatically switches on. The compressor automatically switches off, when pressure in the air tank reaches the switch-off pressure.

Hazard of damage to pneumatic components.

Working pressures of the pressure switch set by the manufacturer must not be changed. If the compressor operates at higher working pressure than the switch-on pressure, it means that the air consumption is too high (see Chapter Malfunctions)

Hazard of damage to dryer.

When the dryer is operated at ambient temperature that is higher than maximum operating temperature stated in Technical data, the dryer can be damaged.

Prior to operation it is necessary to connect the housing using output transition piece for compressed air via internal thread G 1/4". Disconnect the connected device only after the reduction of pressure in the pneumatic system to zero!

Required drying performance can only be achieved when the defined operating conditions are followed.

Drying performance decreases and the achieved dew point drops if the dryer is operated at any pressure below the minimum working pressure.

During longer operation of the compressor, the temperature in the surrounding of the compressor increases above 40 °C and the cooling fan switches on automatically. After cooling the space below 32 °C, the fan switches off.

14. SWITCHING ON THE COMPRESSOR

(Fig. 3) The compressor is activated using the main switch on the cabinet (28). The compressor will start and compressed air will fill the air tank. As the compressed air is used, the pressure in the air tank drops to a pre-set level, the compressor switches on and compressed air is supplied to the tank. Once the shut-off pressure is reached the compressor automatically shuts down. The compressor comes back on once the pressure in the air tank decreases to the switching-on pressure.

The switching pressures are checked using the pressure gauge (32). A tolerance of $\pm 10\%$ is acceptable. The pressure in the air tank may not exceed the permitted operating pressure.

Compressor (without dryer) – when first started and put into operation, the compressor fills the air tank to the switch-off pressure and the pressure switch turns off the compressor. Then the compressor operates in automatic mode, switches on and off depending on compressed air consumption.

Compressor with dryer – the compressor operates in the same manner as above, only the compressed air passes through the cooler that removes moisture from the compressed air.

Risk of damage to pneumatic parts.

The pressure switch (1) was set up by the manufacturer and any other changes to its settings may be performed only by a qualified technician trained by the manufacturer.

15. SWITCHING OFF THE COMPRESSOR

(Fig. 3) The main switch (28) and pulling the mains plug is used to switch off the compressor to perform service or for any other reason. This disconnects the compressor from the power

Vent the pressure in the air tank to zero by opening the drain valve.

PRODUCT MAINTENANCE**16. PRODUCT MAINTENANCE**

The operator should carry out device checks regularly in the intervals defined by applicable regulations. Test results must be recorded.

The device has been designed and manufactured to keep its maintenance to a minimum. For correct and reliable operation of the compressor perform the following operations:

Risk of servicing by persons without required qualification.

Repair works beyond standard maintenance (see chapter.16.1) may only be performed by a qualified technician (organization authorized by the manufacturer) or manufacturer's customer service.

Standard maintenance works (see chapter 16.1) may only be carried out by a trained operator.

Only use spare parts and accessories approved by the manufacturer.

Risk of injury and damage to the device.

Before any maintenance work, it is necessary to:

- check, if it is possible to disconnect compressor from the supplied equipment, so that there is no risk of damage to a person using the given appliance or any other material damages
- switch off the compressor
- disconnect it from the mains (pull out the mains plug)
- vent the compressed air from the air tank

Risk of injury when venting the compressed air.

Protect eyesight – wear goggles when venting the compressed air from the compressor pneumatic system (air tank).

Risk of burn.

Pump components (head, cylinder, pressure hose) have high temperature during and shortly after compressor operation – do not touch these components!

Let the device cool before any product maintenance, service or connection/disconnection of pressurized air!

The removed grounding conductor during service must be connected back to the original position after completing the service.

16.1. Maintenance intervals

Time interval ^{c)}	50 Hz	60 Hz	Once a day	Once a week	Once a year	Once every 2 years	2000 h 4000 h 6000 h 8000 h 10000 h 12000 h						Chapter	Set of spare parts	Performed by
							1600 h	3200 h	4800 h	6400 h	8000 h	9600 h			
Product operational checks	X											16.2	-	user	
Drain the condensate from the air tank - At high RH ^{b)}	X											16.5	-	user	
- Pour the condensate out of the vessel ^{d)}												16.5	-	user	
Drain the condensate from the vessel - At normal RH ^{b)}			X									16.3	-	Qualified professional	
Leak check connections and inspect device					X							16.4	-	Qualified professional	
Inspection of electrical connections					X							16.9	-	Qualified professional	
Cooler and fan check					X							16.8	025200304-000	Qualified professional	
Filter element replacement in the dryer filter ^{d)}					X							16.6	-	Qualified professional	
Safety valve check					X							16.7	025200126-000	Qualified professional	
Replacing the inlet filter and pre-filter on the DK50 PLUS air pump ^{a)}					X			X				16.7	025200139-000	Qualified professional	
Replacing the inlet filter and pre-filter on the DK50 2V air pump ^{a)}					X			X				16.7	025200150-000	Qualified professional	

a) data is presented in hours; if not possible, then the data is considered in years

b) only for compressors without dryers

c) For 60 Hz compressor models, reduce the time interval in hours by 20%

(2000h./1600h., 4000h./3200h., 6000h./4800h., 8000h./6400h., 10000h./8000h., 12000h./9600h.)

d) only for compressors with dryers

16.2. Check of product operation

- Check pump condition – the pumps must operate evenly without excessive vibrations or noise. If there is any negative result, find the cause of the given condition or call service engineers.
- Check fans operation (visually) – the fans must be operating when the fans are running. If there is a negative result, look for the cause of the given condition or call service engineers.
- Check if the power cord and pneumatic hoses are not damaged. Replace the damaged parts or call service engineers.
- Check ambient temperature – ambient temperature must be below temperature limit (40°C). Improve cooling in the room if the temperature is higher.
- Compressor with air dryer - open the plug on the condensate collecting bottle and drain the condensate,

16.3. Check of pneumatic connections leakage and device inspection

Leakage check

- Check pneumatic connections of the compressor when the compressor is operating – pressurizing.
- Use a leakage analyser or soapy water to check all connections for leaks. Tighten or release the connection where leaks are found.

Device inspection

- Check condition of the compressor pump – even operation, appropriate noise level.
- Check fans operation – fans must be running in the defined cycles of the compressor operation.
- Check function of the temperature switch (B2) – heat the temperature switch to temperature above 40°C (e.g. using a heat gun – Be careful, do not operate with high temperature on plastic components in the surrounding, it can cause deformation of plastic materials). Fan EV1 (and EV2 for a compressor with dryer) starts up once the temperature reaches 40°C – the compressor must be under voltage.
- Check filters condition – filters must be undamaged and sufficiently clean.
- Check condition of the pump itself, check if there is no dirt inside the crankcase or clearance in the crankshaft.
- Replace any defective parts if there are any failures.

16.4. Check of electrical connections

Risk of electric shock.

Check electrical connections on the device disconnected from the mains.

- Check mechanical function of the main switch.
- Check if the power cable and conductors are not damaged.
- Visually check if cables are connected to the terminal box.
- Check all screw connections of the green-yellow PE grounding conductor.

16.5. Condensate drain

Risk of slipping on a wet floor in case of the collecting vessel overflow.

Switch off the compressor every time before emptying the vessel!

Carefully and slowly open the drain valve. Quickly opening the drain valve is accompanied by a loud sound and uncontrollable stream of accumulated condensate.

Compressors without air dryers

(Fig. 1 and Fig. 5) During regular operation is condensate automatically drained through automatic drain (autodrain) which is catched into a condensate collection vessel. Pull the vessel out of the holder and pour out the condensate.

Compressors with air dryer

- Condensate from compressors with air dryers is automatically drained into a condensate collection vessel.
- The bottle must be drained regularly.

Fig. 5: Check of condensate collection vessel

Monitor the level in the vessel using the 1 l or 2 l markings (depending on the volume of the vessel), and empty at least once a day. The vessel may overflow if the condensate is not drained out of the vessel at the defined interval.

If necessary, it is possible to drain the air tank using the drain valve (1) from the lower side of the air tank (Fig. 6).

Switch off the compressor at the mains and reduce the air pressure in the device to a maximum of 1 bar, e.g. by venting air through the connected equipment.

Use a hose to drain with the free end placed in a suitable container and then open the valve to drain condensate from the tank.


Fig. 6: Condensate drain

Prior to the following checks it is necessary to disassemble the cabinet door:

16.6. Check of safety valve

When the compressor is operated for the first time, make sure that the safety valve is working properly.

- Turn the screw (2) on the safety valve (1) several rotations to the left until the safety valve releases air.
- Let the safety valve vent for only a few seconds.
- Turn the screw (2) to the right until it seats, closing the valve.


Fig. 7: Check of safety valve

Risk of pressure increase when safety valve is damaged.

Safety valve must not be used to depressurize the air tank. It could damage the safety valve function. The valve is set to the maximum allowed pressure, it is tested and marked.

Safety valve setting must not be adjusted.

Risk of injury when compressed air is vented.

Protect your eyesight when checking the safety valve – wear protective goggles.

16.7. Inlet filter and pre-filter replacement

The inlet filter (1) is located inside the compressor crankcase cover.

Inlet filter replacement:

- Pull out the rubber plug by hand (2).
- Remove the used filter (1).
- Insert a new filter and replace the rubber plug.

Pre-filter replacement:

- Pull out the pre-filter by hand (3).
- Replace the element and re-insert.




Fig. 8: Inlet filter and pre-filter replacement

16.8. Replacement of filter element

- Remove the hose (1) from the quick connector.
- Use a wrench (2) to release the filter vessel (3) and remove.
- Pull down on the filter element (4) to remove.
- Insert a new filter element.
- Re-install the filter vessel.
- Gently tighten the filter vessel with the wrench.
- Re-install the hose on the quick connector.

Fig. 9: Replacement of filter element

16.9. Check of fan and cooler

The equipment, in particular the compressor fan (14), cooler fan, and the cooler (10), must be kept clean to ensure efficient drying. Remove dust from surface of the fan or cooling fins by vacuuming or blowing down with compressed air.(Fig. 1, Fig. 2)

16.10. Cleaning of the exterior surfaces of the product

Use neutral agents for cleaning of the external surfaces of the product.

The use of aggressive cleaning agents and comprising alcohol and chlorides may lead to the damage of the surface and the discoloration of the product.

17. LONG-TERM SHUTDOWN

When there is a predisposition that the compressor shall not be used for a prolonged time period, release air pressure in the pressure tank. Switch off the compressor by switch (28) (Fig. 3) and disconnect the device from the mains

18. DISPOSAL OF DEVICE

- Disconnect the device from the mains.
- Release air pressure in the pressure tank by opening the drain valve.
- Dispose of the device following all applicable environmental regulations.
- Entrust a specialised company to sort and dispose of waste.
- Worn out components have no negative environmental impact.

TROUBLESHOOTING**Risk of electric shock.****Before any of the following operations on the device, disconnect the device from the mains (pull out the mains plug).****Risk of injury during work with pneumatic components under pressure.****Before any of the following operations on the device, it is necessary to decrease pressure in the air tank and in the pneumatic system to zero.****Any operations concerning the troubleshooting may be performed only by a qualified technician.**

Malfunction	Possible cause	Solution
Compressor does not switch on	No voltage in the pressure switch	Check voltage in socket
		Check circuit breaker switch – switch to position switched-on „I“
		Loose conductor from terminal - repair
		Check electrical cord – replace defective cord
	Motor winding failure, damaged thermal protection	Replace motor or windings
	Capacitor failure	Replace capacitor
Compressor often switches on	Seized piston or other rotating part	Replace damaged components
	Pressure switch does not switch	Check function of pressure switch
	Air leakage in pneumatic system	Check pneumatic system – seal loose connections
	Non-return valve (SV) leakage	Clean non-return valve, replace seals, replace non-return valve
Low pressure in the air tank (compressor running constantly)	Large volume of condensed liquid in pressure tank	Drain condensed liquid
	Low compressor capacity	Check time of filling the air tank
	High air consumption of supplied equipment	Decrease air consumption Use compressor with higher capacity
	Leakage in pneumatic system	Check pneumatic system – seal loose connections
	Low capacity of the pump	Clean / replace the pump
Prolonged operation of the compressor	Pump failure	Clean / replace the pump
	Dryer failure	Replace dryer
	Air leakage in pneumatic system	Check pneumatic system – seal loose connection
	Worn piston ring	Replace worn piston ring
	Dirty inlet filter	Replace dirty filter with new filter
	Incorrect function of solenoid valve	Repair or replace fan or coil

Compressor is noisy (knocking, metal noises)	Damaged piston bearing, piston rods, motor bearing	Replace damaged bearing
	Loose (cracked) dampening element (spring)	Replace damaged spring
Dryer does not dry (condensed water in the air) *	Defective cooler fan	Replace fan Check power source
	Damaged dryer	Replace dryer
	Defective automatic condensate drain	Clean / replace

)* Clean inside surfaces of the air tank thoroughly and remove all condensed liquid after dryer failure.

To protect the connected equipment from any damage, it is necessary to check humidity of the released air from the air tank (see chapter Technical data)

19. REPAIR SERVICE

Guarantee and post-guarantee repairs are provided by the manufacturer or organizations and technicians approved by the manufacturer.

Warning

The manufacturer reserves the right to make changes on the device, which will not significantly affect properties of the device.

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ANNEX / ANHANG / PRÍLOHA**20. INSTALLATION RECORD**

1. Product: (model) DK50 PLUS MOBILE DK50 PLUS/M MOBILE DK50 2V MOBILE DK50 2V/M MOBILE		2. Serial number:	
3.1. User's name:			
3.2. Address of installation:			
4. Equipment connected to the compressor:			
5. Installation / Commissioning:		6. Contents of operator training:	
Product completeness check **	Y	Description of the product and functions**	Y
	N		N
Documentation completeness check **	Y	Product operation: turning on/off, controls, control procedures, data on the display panel, alarms, operation in alarm conditions**	Y
	N		N
Installation/connection to equipment **	Y	Product maintenance: maintenance intervals, maintenance procedure, service intervals, operating activities**	Y
	N		N
Functional test **	Y	Safety measures, warnings – their meaning and compliance **	Y
	N		N
Notes:			
7. Operator instructed on safety measures, operations and maintenance:			
Name:		Signature:	
Name:		Signature:	
Name:		Signature:	
8. Installation and instruction performed by: First name/Last name		Signature:	
Company:		Address:	
Phone:			
Email:		Date:	
9. Distributor:			
Company:		Address:	
Contact person:			
Phone:		Email:	

** mark with an "X" in points 5 and 6 (Y - yes /N - no). Enter any observations from points 5 and 6 into the "Notes" section

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