DO 3

- EN User manual
- DE Benutzerhandbuch
- FR Manuel de l'utilisateur
- RU Руководство пользователя
- PL Instrukcja obsługi
- SK Návod na použitie
- CS Návod k použití

DENTAL SUCTION SYSTEM
semi-dry
ZAHNÄRZTLICHES ABSAUGSYSTEM
Halbtrocken
SYSTÈME D'ASPIRATION DENTAIRE
semi-sec
СТОМАТОЛОГИЧЕСКАЯ ВАКУУМНАЯ СИСТЕМА
Полусухая
STOMATOLOGICZNY SYSTEM SSANIA
częściowo na sucho
DENTÁLNY ODSÁVACÍ SYSTÉM
Polosuchý
DENTÁLNÍ ODSÁVACÍ SYSTÉM
Polosuchý

DO 3

EKOM spol. s r. o.
Priemyselná 5031/18
SK-921 01 Piešťany
Slovak Republic
tel.: +421 33 7967255
fax: +421 33 7967223

www.ekom.sk
email: ekom@ekom.sk

DATE OF LAST REVISION
DATUM DER LETZTEN ÜBERARBEITUNG
DATE DE LA DERNIÈRE RÉVISION
ДАТА ПОСЛЕДНЕЙ РЕДАКЦИИ
DATA OSTATNIEJ AKTUALIZACJI
DÁTUM POSLEDNEJ REVÍZIE
DATUM POSLEDNÍ REVIZE

04/2024

NP-DO3-16_04-2024
112000438-000

EN	<u>CONTENTS</u>	<u>5</u>
DE	<u>INHALT.....</u>	<u>49</u>
FR	<u>TABLES DES MATIÈRES</u>	<u>94</u>
RU	<u>ОГЛАВЛЕНИЕ</u>	<u>140</u>
PL	<u>SPIS TREŚCI.....</u>	<u>186</u>
SK	<u>OBSAH.....</u>	<u>231</u>
CS	<u>OBSAH.....</u>	<u>275</u>

CONTENTS

IMPORTANT INFORMATION	6
1. CONFORMITY WITH THE REQUIREMENTS OF EUROPEAN UNION DIRECTIVES ...	6
2. INTENDED USE	6
3. CONTRAINDICATIONS AND SIDE-EFFECTS.....	6
4. WARNINGS AND SYMBOLS	6
5. WARNINGS.....	7
6. STORAGE AND TRANSPORT CONDITIONS.....	9
TECHNICAL DATA	10
7. TECHNICAL DATA.....	10
PRODUCT DESCRIPTION	14
8. PRODUCT DESCRIPTION.....	14
9. PRODUCT FUNCTION.....	18
INSTALLATION.....	21
10. CONDITIONS FOR USE	21
11. POSITIONING OF THE DENTAL SUCTION SYSTEM	21
12. INSTALLATION OF THE PRODUCT	24
13. CONNECTION DIAGRAMS.....	32
14. COMMISSIONING	37
OPERATION.....	38
15. SWITCHING THE PRODUCT ON	38
PRODUCT MAINTENANCE.....	42
16. PRODUCT MAINTENANCE	42
17. MAINTENANCE INTERVALS	43
TROUBLESHOOTING	47
18. REPAIR SERVICE.....	48
19. LONG-TERM SHUTDOWN	48
20. DISPOSAL OF DEVICE.....	48
ANNEX.....	319
21. INSTALLATION RECORD.....	319

IMPORTANT INFORMATION**1. CONFORMITY WITH THE REQUIREMENTS OF EUROPEAN UNION DIRECTIVES**

This product conforms to the requirements of the Regulation (EU) on medical devices (MDR 2017/745) and is safe for the intended use if all safety instructions are followed.

2. INTENDED USE

Dental suction system DO 3 is a source of vacuum for dental units. The device can be used with dental units equipped with a suction block without a water-air separator.

Risk of injury or equipment damage.

The device must not be used to suction aggressive, flammable and explosive mixtures.

Any other use of the product beyond the intended use is considered as incorrect use. The manufacturer is not responsible for any damages or injuries as a result of incorrect use or disobedience to instructions stated in this User manual. All risks shall be solely borne by the user/operator.

3. CONTRAINDICATIONS AND SIDE-EFFECTS

There are no contraindications or side-effects known.

4. WARNINGS AND SYMBOLS

The following symbols are used in the user manual, device and its packaging to denote important details and information:

General caution

General warning

Warning – risk of electric shock

Warning – hot surface

Danger of biological hazard

Warning - controlled automatically

Refer to instruction manual

CE-marking

Medical device

Serial number

IMPORTANT INFORMATION

	Article number
	Unique Device Identifier
	Protecting earthing
	Package handling label – fragile
	Package handling label – this side up
	Package handling label – keep dry
	Package handling label – temperature limits
	Package handling label – limited stacking
	Package label – recyclable material
	Manufacturer

5. WARNINGS

The product is designed and manufactured to be safe for the user and the surrounding environment when used in the defined manner. Keep the following warnings in mind.

5.1. General warnings

READ THE USER MANUAL CAREFULLY BEFORE USING THE DEVICE AND KEEP IT FOR FUTURE USE!

- The user manual supports the correct installation, operation and maintenance of the product. Careful review of this manual will provide the information necessary for the correct operation of the product for its intended use.
- Original packaging should be retained for the return of the device. Only the original packaging ensures protection of the device during transport. If it is necessary to return the product during the warranty period, the manufacturer is not liable for damages caused by improper packaging.
- This warranty does not cover damages originating from the use of accessories other than those specified or recommended by the manufacturer.
- The manufacturer only guarantees the safety, reliability and function of the device if:
 - installation, new settings, changes, expansion, and repairs are performed by the manufacturer or an organisation authorised by the manufacturer,
 - the product is used pursuant to the user manual.
- The user manual corresponds to the configuration of the product and its compliance with the applicable safety and technical standards at the time of printing. The manufacturer reserves all rights for the protection of its configuration, methods and names.

- Translation of the user manual is performed in accordance with the best available knowledge. The Slovak version is to be used in the event of any uncertainties.
- This user manual is the original instructions. Translation is performed in accordance with the best available knowledge.

5.2. General safety warnings

The manufacturer designed and manufactured the product to minimise all risks when used correctly for the intended use. The manufacturer considers it its obligation to lay down the following general safety precautions.

- Use and operation of the product must comply with all laws and local regulations valid in the place of use. The operator and user are responsible for following all the appropriate regulations in the interests of performing work safely.
- Only the use of original parts guarantees the safety of operating personnel and the flawless operation of the product itself. Only the accessories and parts mentioned in the technical documentation or expressly approved by the manufacturer should be used.
- The operator must ensure that the device is functioning correctly and safely before every use.
- The user must be familiar with the operation of the device.
- Do not use the product in environments with a risk of explosion.
- The user must inform the supplier immediately if any problem directly related to the operation of the device occurs.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

5.3. Safety warnings on protection from electric current

- The device must only be connected to a properly installed, earthed socket.
- Before the product is plugged in, ensure that the mains voltage and frequency stated on the product are in compliance with the values of the mains.
- Prior to putting the device into operation check for any damage to the connected pneumatic lines and electrical wiring. Replace damaged pneumatic lines and electrical wirings immediately.
- Immediately disconnect the product from the mains (remove power cord from the socket) in hazardous situations or in the case of a technical malfunction.
- During all repairs and maintenance, ensure that:
 - the mains plug is removed from the power socket,
 - piping, water-air separator and storage tank drained.
- The product shall only be installed by a qualified technician.

IMPORTANT INFORMATION

6. STORAGE AND TRANSPORT CONDITIONS

The manufacturer provides sufficient protection of the appliance when it leaves the factory. Therefore the product is protected against damage during storage and transport.

Risk of liquid spillage.

During transport, the separation container and the storage container have to be empty. The inlets of the amalgam separator have to be closed off with plugs.

Caution! For transport, always use the original packaging and, transport the appliance in the upright position, using secure fixings.

Keep dry. Protect appliance against humidity and extreme temperatures during storage and transport. Appliance in the original packaging must be stored in warm, dry and dust-free areas.

Keep packaging material if possible. If not, please dispose of the packaging material in an environmentally friendly way. Cardboard can be recycled.

Storing or shipping the equipment in any conditions other than those specified below is prohibited.

6.1. Ambient conditions for storage and transport

Products can be stored in rooms and means of transport that are free from any traces of volatile chemical substances under the following climatic conditions:

Temperature	from -25°C o $+55^{\circ}\text{C}$, in 24h to $+70^{\circ}\text{C}$
Relative humidity	10% to 90 %, (without condensation/ non-condensing)

TECHNICAL DATA
7. TECHNICAL DATA


Products are designed to operate in dry, ventilated and indoor dust-free rooms with the following climatic conditions:

Temperature	from +5°C to +40°C
Relative humidity	max. 70%

Type		DO 3		
		BASIC	STANDARD	ADVANCED
Performance measured on the coupling of the large cannula	l/min	300	600	900
Rated voltage	V, Hz	230	230	230
Frequency		50/60	50/60	50/60
Number of connected dental units		1	1-2	3-4
The coefficient of use with simultaneous operation of dental units	%	100	100	75
Maximum vacuum	mbar	-130	-160	-200
Rated current	A	3/4	3.9/5.4	11/14.5
Noise level	dB(A)	60	62	64
Noise level in the cabinet	dB(A)	54	55	-
Operating mode		intermittent, S3	continuous, S1	continuous, S1
Physical data – without/with cabinet				
Height	mm	520 / 545	520 / 545	590
Depth	mm	420 / 445	420 / 445	540
Width	mm	360 / 405	360 / 405	360
Weight	kg	21 / 32	25 / 35	41
Physical data with pre-assembly for amalgam separator – without/with cabinet				
Height	mm	820 / 845	820 / 845	890
Depth	mm	420 / 445	420 / 445	540
Width	mm	360 / 405	360 / 405	360
Weight	kg	25 / 41	29 / 44	46
Requirements for the amalgam separator				
Minimum water flow	l/min	2.5		5
Dental suction system type acc. to ISO 10637		semi-dry		
Classification pursuant to EN 60601-1		Class I		
Classification pursuant to MDR 2017/745		I		

TECHNICAL DATA

7.1. Characteristic curves of the suction pump

7.2. Electromagnetic compatibility declaration

Medical equipment needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided below.

Guidance and manufacturer's declaration - electromagnetic emissions		
Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances		
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The equipment uses RF energy only for its internal functions. Therefore, the RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The equipment is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	The equipment is not likely to cause any flicker, as the current flow is approx. constant after the start up.	

Guidance and manufacturer's declaration - electromagnetic immunity			
Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances			
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV 100 kHz repetition frequency Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV L-N ±2 kV L-PE; N-PE Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruption, and voltage variations on power supply input lines IEC 60601-4-11	$U_T=0\%$, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°) $U_T=0\%$, 1 cycle $U_T=70\%$ 25/30 cycles (at 0°) $U_T=0\%$, 250/300 cycles	$U_T=>95\%$, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°) $U_T=>95\%$, 1 cycle $U_T=70\%$ (30% dipp U_T), 25(50Hz)/30(60Hz) cycles (at 0°) $U_T=>95\%$, 250(50Hz)/300(60Hz) cycle	Mains power quality should be that of a typical commercial or hospital environment. The device stops and restarts automatically at each dip. In this case does not occur unacceptable pressure drop.
NOTE: U_T is the A.C. mains voltage prior to application of the test level.			

TECHNICAL DATA

Guidance and manufacturer's declaration - electromagnetic immunity			
Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances			
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1,2\sqrt{P}$ $d=1,2\sqrt{P}$, 80 MHz to 800 MHz $d=2,3\sqrt{P}$, 800 MHz to 2,7 GHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz	3 V/m	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Proximity fields from RF wireless communications equipment IEC 61000-4-3	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment.			
^b Over the frequency range 150 kHz to 80 kHz, field strengths should be less than 3 V/m.			

PRODUCT DESCRIPTION**8. PRODUCT DESCRIPTION**

The DO 3 is a semi-dry dental suction system with a medium flow. Included is a water-air separator, some variants are ready for additional installation of the amalgam separator (the amalgam separator is not included). In order to reduce the noise levels, it can be equipped with a cabinet.


It is suitable for dental units equipped with a suction block without a water-air separator. Suction system models with pre-assembly for the amalgam separator are designed for cases where the amalgam separation is required. The design of the suction system makes it suitable to be placed in the technical background of the dental work site. If it is equipped with a soundproof cabinet, it may be placed in the dental work site.

8.1. Variants


Dental suction systems come in the following variants according to performance and range of accessories:

Variant	Version	Description
DO 3	BASIC	Suction system with a performance of 600 l/min., designed for 1 dental unit. It is placed in the soundproof cabinet. It operates in an intermittent mode - the aspirated liquid has to be emptied during downtime when the suction is turned off. The suction system cannot be equipped with an amalgam separator.
DO 3 Z	BASIC	DO 3 BASIC variant without cabinet.
DO 3 A	BASIC	DO 3 BASIC with pre-assembly for amalgam separator.
DO 3 ZA	BASIC	DO 3 BASIC with pre-assembly for amalgam separator, without cabinet.
DO 3	STANDARD	Suction system with a performance of 600 l/min., designed for 1-2 dental units, suitable for continuous operation without downtime. It is placed in the soundproof cabinet. It can be equipped with an amalgam separator.
DO 3 Z	STANDARD	DO 3 STANDARD variant without cabinet.
DO 3 A	STANDARD	DO 3 STANDARD with pre-assembly for amalgam separator.
DO 3 ZA	STANDARD	DO 3 STANDARD with pre-assembly for amalgam separator, without cabinet.
DO 3 Z	ADVANCED	Suction system with a performance of 900 l/min., designed for 3-4 dental units, suitable for continuous operation without downtime. It can be equipped with an amalgam separator. The suction system is not equipped with a cabinet.
DO 3 ZA	ADVANCED	DO 3 ADVANCED with pre-assembly for amalgam separator, without cabinet.


PRODUCT DESCRIPTION


DO 3 Z BASIC
DO 3 Z STANDARD
DO 3 Z ADVANCED

DO 3 BASIC
DO 3 STANDARD

DO 3 ZA BASIC
DO 3 ZA STANDARD
DO 3 ZA ADVANCED

DO 3 A BASIC
DO 3 A STANDARD

Meaning of the letters in the designation:

Z - without cabinet

A - with pre-assembly for amalgam separator







8.2. Accessories

Accessories are included in the standard delivery.



Type	Use	Article number	
Exhaust hose DN40 mm 1.5 m	BASIC STANDARD ADVANCED	062001393-000	
Hose clamp 25-40 4 pcs	BASIC STANDARD ADVANCED	049000263-000	
Nozzle DN32/30 1 pc	BASIC STANDARD	024003049-000	
Nozzle DN40/40 1 pc	ADVANCED	062001489-000	
Reducer DN32/25 1 pc	BASIC STANDARD ADVANCED	604041206-000	
Separator filter (Dust filter)	BASIC STANDARD ADVANCED	025200454-000	

8.3. Accessories

Accessories that are not included in the standard order must be ordered separately.

Type	Use	Article number	
Exhaust hose 3 metres, 2 clamps DN40 mm	BASIC STANDARD ADVANCED	604041200-000	
Cabinet	BASIC STANDARD	603013644-000	
Cover for the amalgam separator	BASIC STANDARD	603013991-000	
Amalgam separator frame	BASIC STANDARD	604022817-000	
	ADVANCED	604022805-000	
Amalgam separator connections	BASIC STANDARD ADVANCED	604041252-000	
Amalgam Separator Kit	STANDARD	604041239-000	

PRODUCT DESCRIPTION

Amalgam Separator Kit	ADVANCED	604041240-000	
Wall mounting kit	BASIC STANDARD	604041202-000	
Vacuum valve DN32	BASIC STANDARD	604022820-000	
Vacuum valve DN40	ADVANCED	604022821-000	
Silencer with bacterial HEPA air filter	BASIC STANDARD	604032022-000	
	ADVANCED	604032103-000	
Bacterial air filter element (class E12 per EN1822)	BASIC STANDARD	025200426-000	
	ADVANCED	025200468-000	
Extension	BASIC STANDARD ADVANCED	604032181-000	

9. PRODUCT FUNCTION

(Fig. 1): The source of vacuum in the suction pipe is a suction pump (16). The air aspirated from the dental units containing liquids and solids is transported from the main suction pipe through a flexible hose (18), particles larger than 2 mm are trapped on the screen of the inlet filter (11). The liquids and solids are separated from the air in a water-air separator (10) and pumped by a drain pump (12) into the storage container (7). When the amalgam separator is retrofitted, the accumulated liquid containing solid particles flows freely through the amalgam separator (14), which traps solid particles containing amalgam. The waste liquid flows through the hose (19) into the sewer. Air from the water-air separator (10) passes through a suction pump (16) and enters the exhaust pipe connected via a flexible hose at the coupling (15). The pipe is connected with the exterior through an opening in the cabinet.

The mains cable is led out of the extraction system. The dental unit cord is connected to the corresponding terminal board on the controller (13). The operation of the suction pump is then controlled by closing the contacts when lifting the fitting of the suction cannula of the dental unit. The drain pump is switched based on the signal from the level sensor (8). If the level rises up to the upper level then, based on the signal from the level sensor (9), the suction pump (16) is switched off as a result.

The soundproof cabinet (4) reduces the amount of noise that would otherwise spread into the surrounding area. Forced air ventilation is provided by a fan (17) under the suction pump (16). The cooling air is supplied via a slot between the cabinet of the suction system and the floor; it is then expelled through openings in the top part of the cabinet.

When retrofitting the amalgam separator, follow the Amalgam Separator Installation Instructions.

PRODUCT DESCRIPTION

Fig. 1: DO 3



Fig. 2: Cabinet

Descriptions for figures 1 - 2

- | | | | |
|----|---|----|--|
| 1 | Power on (green)** | 12 | Drain pump ** |
| 2 | Warning signal (yellow)** | 13 | Controller |
| 3 | Display | 14 | Amalgam separator location * |
| 4 | Soundproof cabinet | 15 | Exhaust hose connection point |
| 5 | Mains switch | 16 | Suction pump |
| 6 | Cover for amalgam separator * | 17 | Fan |
| 7 | Storage container * ** | 18 | Hose for connection to the main suction pipe |
| 8 | Lower level sensor** | 19 | Hose for connection to the sewer pipe |
| 9 | Upper level sensor | | |
| 10 | Water-air separator | | |
| 11 | Housing for the screen of the inlet filter ** | | |
- *) applicable only to the model with an amalgam separator
 **) does not apply to BASIC models

INSTALLATION

INSTALLATION

10. CONDITIONS FOR USE

- The appliance must be installed and operated in a dry, well ventilated and dust-free area where ambient temperature is within the range of +5°C to +40°C and relative air humidity does not exceed 70%. The appliance must be installed so that it is accessible at all times to operate and maintain. Please ensure the appliance label is accessible also.
- The device may be installed on the same level as dental units, in a side room or on a lower level.
- The device must stand on a flat, sufficiently stable base (be aware of the compressor weight, see Chapter 7. Technical data).
- Installation on the wall is only suitable for building materials with sufficient strength (brick, concrete). It is important to select a proper anchoring type during installation (pay attention to the weight of the device, see Chapter 7. Technical data). The manufacturer is not liable for risks and damages caused by improper anchoring of the product onto the wall.
- Suction system must not be exposed to external environments. They must not be in operation in moist or wet environments. Furthermore, it is forbidden to use them in premises, where there is the presence of explosive gases or flammable liquids.
- Any other use or use outside this framework is not considered as intended use. The manufacturer is not responsible for any damages arising from such use. The operator/user assumes all risks.

11. POSITIONING OF THE DENTAL SUCTION SYSTEM

11.1. Unpacking the product

- Unwrap the suction system from the packaging.
 - Remove the cover
 - Remove the fixation elements
 - Remove the casing and packing material
 - Remove the product

Fig. 3: Unpacking the product without the cabinet





Fig. 4: Unpacking the product with the cabinet

11.2. Handling



- Position the suction system at the operation site (Fig. 5)

Fig. 5: A: Handling without a cabinet





INSTALLATION


B: Handling with a cabinet

disconnect all earthing cables

connect all earthing cables

12. INSTALLATION OF THE PRODUCT

Check all piping connections for leaks after installation. Inspect the functionality of the electric controls from each dental unit and verify the proper functioning of the product as described in Chapter 9.

Risk of incorrect installation.

Only a qualified professional can install and commission the suction system. Their duty is to train operating personnel on the use and maintenance of the equipment. The installation and training of the operating staff have to be confirmed by the installer's signature on the certificate of installation (see Annex).

Danger of burns or fire.

Some parts of the suction system may become hot during operation and reach hazardous temperatures (dangerous to the touch of personnel or materials).

All used pipes or hoses have to be vacuum resistant and resistant to all chemicals which are normally used in dental practice (e.g. sewer pipes made from PP, PVC-C, PVC-U, PE-HD).

The connection of the dental suction system DO 3 to the installation pipes must be made with flexible hoses that should be as short as possible.

The hose couplings have to be secured with hose clamps.

When replacing the amalgam separator, the clock has to be reset by pressing the separator button located on the display - hold the button pressed for 10 seconds (see chapter OPERATION).

Risk of overheating of the suction system.

Never cover the inlets for the cooling air on the cabinet (around the bottom of the cabinet) and the hot air outlets on the top of the cabinet.

If the device is placed on a floor with a soft surface (e.g. carpet) it is necessary to create a gap for the passage of cooling air between the base-floor and the cabinet-floor e.g. by placing the feet on hard pads.

Risk of damage to pneumatic components.

Suction tubings must not be broken.


You may notice a "new product" odour when you first place the product into service (for a short period of time). This odour is temporary and does not impede the normal use of the product. Ensure the space is properly ventilated after installation.

INSTALLATION

12.1. Installation without the cabinet


Fig. 6: Hose installation

1. Exhaust air hose
2. Hose for connection to the sewer pipe
3. Hose for connection to the main suction pipe

When installing the device on the floor, use the entire rear of the base or openings at the bottom of the base to connect the hose to the main suction line (a), the sewer line (b), and the cable routing.

Fig. 7: Installing on the floor

12.2. Installation with the cabinet

The cabinet can only be used with the models BASIC and STANDARD installed on the floor.

- Prior to the installation, it is necessary to take out the housing with the screen, separate the lid (c) by pulling upwards, and pull the cabinet upwards.

- The suction system is switched on by switching on the switch when using one of the suction hoses with cannula. The dental unit control signal to the suction system is carried by a two-conductor cable. The terminal board for connecting the control cable is located in the controller electrical box (13) (Fig. 1). Connect as shown in the diagram (Fig. 13).
- After connecting the pipes and connecting the cables from above, put the cabinet back. Insert the exhaust air hose via the opening in the back of the cabinet (d) and connect it at the connection point (e).
- Secure the hose with a hose clamp.
- Connect the grounding cable between the suction system and between the cabinet and the cover.
- Press down gently to place the cover onto the cabinet and place the housing with the inlet filter screen into the opening in the cabinet.

INSTALLATION

Fig. 8: Installation with the cabinet


Fig. 9: Connecting the control cable




Diagram for connecting the dental unit control signal:

- Connection via a switch
- 24 V DC wiring

INSTALLATION


12.3. Installation into the wall

Risk of injury due to improper installation on the wall.

It is important to select a proper anchoring type for the material during installation onto the wall (pay attention to the weight of the device, see Chapter 7. Technical data). The manufacturer is not liable for risks and injuries caused by improper anchoring of the product onto the wall.

Fig. 10: Installation into the wall


Piping and connection of the suction system to dental units should be carried out according to the electrical and installation diagrams (see chapter 13) and in accordance with the rules described below.

12.4. Coupling of the main suction pipe

- When multiple dental units are connected, these must be fitted with a local valve.
- Connect the flexible hose (18) to the main suction pipe (Fig. 1). The hose has to have a smooth inner surface with an internal diameter of DN 30 mm (BASIC, STANDARD) or DN 40 mm (ADVANCED). Secure the connection with a hose clamp.
- In order not to lose any suction power, the inner diameter of the main suction pipe has to be at least 32 mm (BASIC, STANDARD) or 40 mm (ADVANCED). The total length of the suction pipe should be as short as possible with a maximum of 15 m (BASIC, STANDARD) or a maximum of 20 m (ADVANCED). The lowest height of the suction pipe must not be lower than 0.4 m below the suction system level.
- Avoid bends at right angles, use bends of 2 x 45° instead.
- We recommend installing the vacuum valve for each dental unit as far as possible in the suction branch (Chapter 13.2).

Fig. 11: Hose connections

1. Reducer 32/25 (connection to the sewerage system)
2. Mouthpiece (connection to suction pipe)
18. Hose (DN30 mm) for connection to the suction pipe (BASIC, STANDARD) or hose (DN40 mm) for connection to the suction pipe (ADVANCED)
19. Hose (DN25mm) for connection to the sewer pipe

12.5. Coupling for the exhaust pipe


Risk of damage due to high temperature.

Use only hoses and pipes made from heat-resistant materials ($\geq 130\text{ }^{\circ}\text{C}$) to create a connection for the exhaust air.

- Connect the exhaust pipe to the connection point (15) (Fig. 1) by using a flexible hose with a smooth inner surface. The diameter of the exhaust air hose connection has to be 40 mm.
- Due to the possibility of the occurrence of an undesirable odour, the exhaust pipe outlet has to be located in the exterior above the level of the roof of the technical background and of the building and away from windows and other air inlets. The outlet has to be designed to prevent particles from entering. If necessary, a bacterial filter with a silencer must be connected in the exhaust pipe (see the installation diagrams).
- The exhaust outlet of the suction system must not be located in the technical room.
- The exhaust outlet has to be designed to prevent liquids, solids and rainwater from entering the system.
- Precautions have to be taken so that the exhaust air does not reach the air intake inlet of the compressor.
- The exhaust pipe has to be routed separately from the connection point of the exhaust pipe to the outlet point outside the building. It should not be connected to any other pipe, e.g. the HVAC pipe.

INSTALLATION

Fig. 12: Examples of a bacterial filter installation in the exhaust pipe

Description to Figure

a	Ventilation grille	
b	Silencer with bacterial filter	article number: 604032022-000
c	Extension	article number: 604032181-000

12.6. Sewer coupling

- Connect the hose for connection to the sewer pipe (19) to the sewer pipe with the reducer (1) (Fig. 11). Provide obstruction-free drainage of the sewage water from the hose.
- The sewer pipe must comply with the relevant local regulations or DIN 1986, parts 1 and 2.
- Provide obstruction-free drainage of the sewage water from every part of the pipe. The sewage pipes including the flexible hose must have a fall of at least 2°. (Chap. 13.2)

12.7. Electrical connection

Fire hazard and risk of electric shock.

The power cord for connecting to the mains must not be broken or damaged.

Risk of electric shock.

Plug the electrical cord into the mains.

The device is equipped with a grounded plug. Ensure full compliance with all local electrical codes. The mains voltage and frequency must comply with the data stated on the device label.

- Keep the socket easily accessible to ensure that in an emergency, the device can be safely disconnected from the mains.
- The connection to the power distribution box must be limited by a current protection device 10A (BASIC, STANDARD) or 16A (ADVANCED).

13. CONNECTION DIAGRAMS


13.1. Electric diagram

Fig. 13: Electric diagram DO 3

1/N/PE~230V, 50/60Hz

~150V, 60Hz

ELECTRICAL OBJECT OF 1st. CAT.

INSTALLATION


Description to electrical diagram

M1	Suction pump	K1	Control electronics
M2	Drain pump	B2	Temperature switch
E1	Fan	S1	Suction switch
F1	Breaker switch	SD	Lower level sensor
X1	Terminal box	SH	Upper level sensor

13.2. Installation diagrams


Fig. 14: Installatio diagrams

A) Installation in the technical room


INSTALLATION

B) Installation in the practice (BASIC, STANDARD)

C) Installation in the basement

INSTALLATION

DO3	BASIC	STANDARD	ADVANCED
 Suction pipe	Min. DN32	Min. DN32	Min. DN40
 Exhaust pipe	Min. DN40	Min. DN40	Min. DN40
 Waste pipe	Min. DN32	Min. DN32	Min. DN40
Ls –maximum length of the suction pipe	15m	15m	20m

Legend for the installation diagrams

- | | |
|-------------------|-------------------|
| a) exhaust outlet | d) terminal strip |
| b) exhaust grid | e) air inlet |
| c) noise muffler | f) vacuum valve |

14. COMMISSIONING

- Check the correct connection of the vacuum air distribution and the sewer connection.
- Check the correct connection of the device to the mains and electrical connection to the dental unit.
- The dental suction system is controlled from the unit; it starts creating a vacuum in the piping after connecting a signal from the unit.

The dental suction system is not equipped with a backup power supply.

OPERATION

The equipment may only be operated by the trained staff!

Risk of electric shock.

In case of emergency, disconnect the device from the mains (pull out the mains plug).

Danger of burns or fire.

Some parts of the suction pump may become so hot during operation that touching them is dangerous to people or material.

The suction pump always activates after the signal is transferred from the dental unit and remains active for approximately 20 seconds after it has been switched off. The cooling fan and drain pump can be operated even when the suction system is switched off.

Risk of overheating the device.

Never cover the inlets for the cooling air on the cabinet (around the bottom of the cabinet) and the hot air outlets on the top of the cabinet.

If the device is placed on a floor with a soft surface (e.g. carpet) it is necessary to create a gap for the passage of cooling air between the base-floor and the cabinet-floor, e.g. by placing the feet on hard pads.

To prevent overheating of the suction system, it is necessary to install a vacuum valve (f -Fig. 14) for the Standard and Advanced versions.

15. SWITCHING THE PRODUCT ON

Switch the mains switch (5) (Fig. 2) to position 1. After switching on, the green power-on signal (1) lights up (Fig. 2). The suction system activates (creates a vacuum in the piping) after the signal is transferred from the dental unit. After the signal is interrupted, the suction system **remains active for approximately 20 seconds**. (factory settings, the interval can be changed – see the service manual)

OPERATION

The status of the device is shown on the touch screen (Fig. 15):





Fig. 15: Display


Meaning of symbols on the display

 1000 Hrs	Total number of operating hours
  	<p>Status of the amalgam separator</p> <ul style="list-style-type: none"> - without amalgam separator - filled 95% or more (yellow) - filled 100% (red) <p>To reset the state of the separator after it has been replaced, press and hold the separator button for 10 seconds.</p>
  	<p>Suction pump</p> <ul style="list-style-type: none"> - switched off - switched on (green) - switched off, liquid in the upper level (red)
  	<p>Drain pump</p> <ul style="list-style-type: none"> - switched off - switched on (green) - not draining (red)
   	<p>Servicing</p> <ul style="list-style-type: none"> - less than 200 hours left before servicing (yellow) - service required (red) - level sensor error
ENG	Language selection
MENU	Open the service menu. Only accessible to service personnel

If the symbols on the display are backlit in red, a warning symbol (2) will light up next to the display (Fig. 2). In this case, the personnel will do the following.

Status of the amalgam separator

The amalgam separator is full. The device remains functional but the amalgam separator may no longer be effective. The personnel are obliged to ensure the separator is collected and recycled.

Suction pump

The suction pump was shut down because of the high liquid level in the water-air separator. If the backlight is red just for a brief moment, it means that the suction system aspirated too much liquid. The suction pump will start again once the aspiration has finished, the device is operational. If the backlight is permanently red, part of the pipe may be blocked. The device is not operational and the personnel have to provide maintenance.

OPERATION

Drain pump

The drain pump does not have the capacity to aspirate the liquid from the water-air separator. This might be caused by a large amount of aspirated liquid, partial blockage of the pipe or by blockage of the pump. If the backlight lights up only temporarily, it might be just an interim effect, the device remains operational. If the red backlight of the symbol is permanent, the personnel are obliged to provide maintenance.

Servicing

The service interval of 5,000 hours has passed since the last service. The personnel have to provide service.

Level sensor error. The operator is obliged to provide service.

PRODUCT MAINTENANCE**16. PRODUCT MAINTENANCE**

The operating entity is obliged to ensure that all tests of the equipment are carried out repeatedly at least once within every 24 months (EN 62353) or in intervals as specified by the applicable national legal regulations. A report must be prepared on the results of the tests (e.g.: according to EN 62353, Annex G), including the measurement methods used.

The device has been designed and manufactured to keep its maintenance to a minimum. For correct and reliable operation perform the following operations.

Risk of servicing by persons without required qualification.

Repair works beyond standard maintenance (see Chapter.17 - Maintenance intervals) may only be performed by a qualified technician (organization authorized by the manufacturer) or manufacturer's customer service. Standard maintenance works (see Chapter 17 - Maintenance intervals) may only be carried out by a trained operator. Only use spare parts and accessories approved by the manufacturer.

Warning – automatic start.

The suction pump always activates once a signal is received from the dental unit and remains active for approximately 20 seconds after it has been switched off. The cooling fan and drain pump can be operated even when the suction system is switched off.

Danger of injury or equipment damage.

Prior to commencing device maintenance, it is necessary to:

- check if it is possible to disconnect the device from the appliance in order to avoid any risk of injury to the person using the appliance or other material damage;
- check if there are no residues of aspirated liquids in the device and drain any residual liquids into the sewer system;
- turn off the device via the mains switch;
- disconnect from the mains (pulling the cord out of the mains socket).

Burn hazard.

Let the device cool off before carrying out service and maintenance!

The removed grounding conductor during service must be connected back to the original position after completing the service

Biohazard warning.

Parts of the device may be contaminated. Before repairing them, it is necessary to have these parts decontaminated by a specialised company.

PRODUCT MAINTENANCE

17. MAINTENANCE INTERVALS

Time interval ^{a)}	2 x day	once a week	once a year		once every 2 years	Chapter	Set of spare parts	Performed by
				5000 h				
Cleaning and disinfection	X					17.3	-	operator
Cleaning the screen of the inlet filter		X				17.4	-	
Cleaning/replacing of the dust filter ^{b)}		X				17.7	025200454-000	
Replacing the bacterial HEPA filter cartridge ^{a)}			X			17.5	025200426-000	qualified technician
Checking connections and visual inspection of the device			X			17.1	-	
Conduct the " Repeated test" as per EN 62353				X		16	-	

^{a)} data is presented in hours; if not possible, then the data is considered in years

^{b)} adapt the cleaning interval to the intensity of use of powdered materials in dental hygiene

17.1. Checking the compressed air connections for leaks and inspection of the equipment

Checking the tightness of the liquid circuit:

- The blower is started and kept running until water enters the liquid circuit. All connections in the circuit are checked.
- There must be no traces of moisture on the joints. Reseal leaky joints – tighten the clamping tapes or replace seals

17.2. Inspection of electrical connections

Inspect the product's electrical connections when the mains are disconnected!

Inspection

- Check the mechanical operation of the main switch.
- Check to ensure the power cable and connected wires are undamaged.
- Check all screw connections of the protective green-yellow PE conductor.

17.3. Cleaning and disinfection

For cleaning and disinfection, use only products approved for dental system systems that do not contain halogens, organic compounds and chlorine.

Disinfection:

- Rinse the cuspidor bowl briefly after each patient completes treatment.
- After each patient completes treatment, flush all suction hoses with water (or disinfectant).
- We recommend using a disinfectant for dental suction systems twice a day (once every 4 hours after being placed into service on a given day).
- The cuspidor bowl must be rinsed twice a day using a disinfectant for dental suction systems.

Cleaning:

Cleaning procedure:

- Prepare 1 litre of cleaning solution in the concentration recommended by the manufacturer into a container.

Caution: do not exceed the prescribed volume, as any larger volume of the cleaning solution drawn in will be automatically drained from the suction system!

- Aspirate the cleaning solution through all cannulae and leave the suction running during the product treatment period (as recommended by the product manufacturer).
- Once the entire volume of cleaning solution has been aspirated, remove the cannulae from the container and do not hang them in the holder so that the cleaning solution can treat the suction system.
- At the end of the treatment period prescribed by the cleaning solution manufacturer, flush the system with 2 liters of clean water.

PRODUCT MAINTENANCE

17.4. Cleaning the housing for the screen of the inlet filter


Biohazard warning.

The screen contains biological waste. Use protective gloves.

- Clean the inlet screen after finishing work in the dental practice and after cleaning and disinfection. All liquids from the piping have to be drained and pumped into the sewer system - the operation symbols of the suction pump and the drain pump must not light up green.
- Switch off the mains switch (1) (Fig. 16).
- Rotate and remove the housing with the screen (2). Unscrew the screen (3) from the housing (4). Empty the content of the screen into the container intended for waste containing amalgam.
- Screw the screen into the housing. Place the housing with the screen back into the filter housing and rotate to tighten.
- Turn on the mains switch.

Recommended cleaning interval of the inlet screen is once per week. The replacement interval can be adjusted in order to secure that the material retained on the screen will not reduce the flow in the suction pipe.


Fig. 16: Cleaning the housing for the screen of the inlet filter

17.5. Bacterial filter element replacement

The bacterial filter is not part of the basic accessories. It may be ordered separately.

Fig. 17: Bacterial filter element replacement

17.6. Cleaning and disinfection of the exterior surfaces of the product

Clean and disinfect the exterior surfaces with neutral cleaning products.

Clean the product with commercial detergents. It is not allowed to clean the product with cleaning products containing an abrasive component, chemical solvents and other aggressive substances.

17.7. Cleaning the dust filter in the separation container

The dust filter is located in the upper lid of the separation container. Its cleaning/replacement interval depends on the intensity of use of powdered materials in dental hygiene: from an interval of once a week for the intensive use of powdered materials in dental hygiene, to an interval of once a year for dental workplaces that do not perform dental hygiene. The filter can be washed in water; in case of mechanical damage, the filter must be replaced with a new one.


- Open the lid of the separation container (Fig. 18)
- Unscrew the white nut (1), remove the grid (2) and pull out the filter (3). (Fig. 18)
- Insert a clean filter (3), fit the grid (2) and screw the nut (1).
- Close the lid of the separation container.

Risk of device overheating. Risk of damage to the device.

When using the device with a clogged filter, the suction power decreases and the device may overheat. When operating with a damaged filter, there is a risk of permanent damage to the device due to dust entering the suction unit.

Fig. 18: Cleaning the dust filter

TROUBLESHOOTING

TROUBLESHOOTING

Risk of electric shock.

Before any of the following operations on the device, disconnect the device from the mains (pull out the mains plug).

Danger of biological hazard and environmental pollution.

To follow the rules of personal hygiene for works with contaminated material. To separate, label, packing and providing for decontamination of contaminated parts by course of national regulations.

Any operations concerning the troubleshooting may be performed only by a qualified technician.

Malfunction	Possible cause	Solution
The suction system is not operational	There is no voltage in the socket	Switch on the circuit breaker, replace the fuse
	Faulty electrical cord	Replace the cord
	Missing control voltage	Check the power cord
		Check the signal from the unit
	The suction system has been switched off by a thermal circuit breaker due to overload, the switch is in position 0	Switch the mains switch ON; if the problem recurs, call for assistance
Damaged motor winding, damaged thermal protection	Replace the motor	
The suction system is not operational, the suction pump symbol is green	Capacitor failure	Replace the capacitor
	The water-air separator is full due to a clogged drain pump or sewer drain	Clean the pump or the clogged drain
The suction system is not operational, the suction pump symbol is red	Drain pump fault	Replace the pump
The suction system is not operational, the drain pump symbol is red	The drain pump/sewer drain is clogged	Clean the pump or the clogged drain
	There is a foreign object in the inlet pipe	Remove the foreign object
The suction system is loud	The suction pump is/has been flooded with liquid	Hand the suction system over to the authorized service provider for repairs - comply with all rules for working with contaminated materials
	The screen of the inlet filter is clogged	Clean the screen of the inlet filter
The suction system does not generate vacuum at all, or generates weak vacuum	Dirty bacterial filter	Change the bacterial filter
	Dirty dust filter in the separation container	Clean/replace the dust filter
	Leaks in the suction pipe	Remove any leaks in the suction pipe
	A foreign object in the inlet/outlet pipe	Remove the foreign object

The key symbol with time indication lights up yellow	Service interval approaching. The time indication contains the number of hours that remain until regular service	Provide service
The key symbol lights up red	The service interval has expired. It is necessary to provide service.	Provide service

18. REPAIR SERVICE

Guarantee and post-guarantee repairs are provided by the manufacturer or organizations and technicians approved by the manufacturer.

Warning.

The manufacturer reserves the right to make changes on the device, which will not significantly affect properties of the device.

19. LONG-TERM SHUTDOWN

If the device will not be used for an extended period of time, it is recommended to clean and disinfect the device in accordance with the procedure in chapter 17.3 and to clean the screen of the inlet filter in accordance with chapter 17.4. Then turn on the mains switch (5) to the "I" position, remove the suction hoses from the dental kit holder and allow air to be drawn in (approx. 15 to 20 min) through the suction hoses so that the suction system dries out completely. Subsequently, after this operation, store the hoses in the holder, turn off the mains switch to the "O" position, and disconnect the device from the mains (pull out the mains plug).

20. DISPOSAL OF DEVICE

Danger of biological hazard and environmental pollution.

The device may be contaminated. Follow personal hygiene rules for handling contaminated materials. When disposing of contaminated components follow local regulations.

- Check that the product is unplugged from the mains.
- Disconnect the product from the pneumatic piping.
- Separate, label, pack and secure decontamination of the contaminated parts in accordance with local regulations.
- Dispose of the product following all locally applicable regulations.
- Decontaminated plastic components of the suction system can be recycled.
- Dispose of electronic circuits and components as you would dispose of electronic waste.
- Other metal components can be disposed of as normal metal scrap.

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21 INSTALLATION RECORD

1. Product: (model) DO 3 Z STANDARD DO 3 ZA STANDARD DO 3 STANDARD DO 3 A STANDARD DO 3 Z ADVANCED DO 3 ZA ADVANCED		2. Serial number:	
3.1. User's name:			
3.2. Address of installation:			
4. Equipment connected to the product:			
5. Installation / Commissioning:		6. Contents of operator training:	
A. Product completeness check **	Y	A. Description of the product and functions**	Y
	N		N
B. Documentation completeness check **	Y	B. Product operation: turning on/off, controls, control procedures, data on the display panel, alarms, operation in alarm conditions**	Y
	N		N
C. Installation/connection to equipment **	Y	C. Product maintenance: maintenance intervals, maintenance procedure, service intervals, operating activities**	Y
	N		N
D. Functional test **	Y	D. Safety measures, warnings – their meaning and compliance **	Y
	N		N
Notes:			
7. Operator instructed on safety measures, operations and maintenance:			
Name :		Signature:	
Name:		Signature:	
Name :		Signature:	
8. Installation and instruction performed by – First name/Last name		Signature:	
Company:		Address:	
Phone:			
Email:		Date:	
9. Distributor:			
Company:		Address:	
Contact person:			
Phone:		Email: :	

** mark with an "X" in points 5 and 6 (Y - yes /N - no). Enter any observations from points 5 and 6 into the "Notes" section

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