

Compressor and wet suction in noise reduction cabinet

- clean oil-free compressed air and vacuum for 1 dental unit
- quiet operation ≤ 52 dB(A)
- more compact than individual compressor and wet suction
- maintenance friendly – toolless lid opening or cabinet removal for easy access to all parts for routine maintenance
- the device can be used for all types of dry or wet suction ready dental units



Technical specification	DBXe	
Nominal voltage, frequency	230, 50/ 60	V, Hz
Compressor efficiency	50/ 60	lit.min ⁻¹ / 5bar
Working pressure	5 - 7	bar
Maximum vacuum	-110	mbar
Suction output at large cannula	300	l.min ⁻¹
Air tank capacity	8,3	l
Mode of operation	S1-100%	
Input power	6,5/ 7,5	A
Noise level	≤ 52	dB(A)
Weight	71	kg
Dimensions (w x l x h)	715x542x680	mm