

Compressor and Suction in one device

- wet or dry suction up to 2 dental chairs
- touch screen with operating parameters
- oil-free compressor with possible output air purity until 0,01 microns
- dry compressed air up to -20°C ADP
- automatic discharging of collected liquids
- air silencer with HEPA filtration
- noise reduction cabinet



Technical specification	DBX	
Nominal voltage, frequency	230/50, 230/60	V, Hz
Compressor efficiency	60 / 135	lit.min ⁻¹ / 6bar
Working pressure	6-8	bar
Air tank capacity	17,5	l
Suction output at large cannula	600	l.min ⁻¹
Maximum vacuum	160	mBar
Mode of operation	S1 / S3	
Input power	7,8 / 14,7	A
Noise level	≤50 / 53	dB(A)
Weight	129 / 146	kg
Dimensions (w x l x h)	600x600x1481	mm