



## Compressor and dry suction in noise reduction cabinet

- oil-free compressed air and vacuum for your dental unit
- quiet operation  $\leq 52$  dB(A)
- smaller footprint than most single compressors
- nonstop operation - advanced cooling allows continual operation
- maintenance friendly - toolless lid opening or cabinet removal for easy access to all parts for routine maintenance
- the device can be used for all types of dry suction ready dental units



Technical specification	DBX-De	
Nominal voltage, frequency	230, 50/ 60	V, Hz
Compressor efficiency	50/ 60	lit.min <sup>-1</sup> / 5bar
Working pressure	5 - 7	bar
Maximum vacuum	-110	mbar
Suction output at large cannula	300	l.min <sup>-1</sup>
Air tank capacity	8,3	l
Mode of operation	S1-100%	
Input power	6,5/ 7,4	A
Noise level	$\leq 52$	dB(A)
Weight	65	kg
Dimensions (w x l x h)	530x550x780	mm