



## Mobile, silent and fully autonomous dental suction in new design

- improved workflow – intuitive manipulation thanks to added top handle
- safer handling – spillage free waste collection vessel
- easier maintenance - hinged door opening, pull-out plate with waste collection vessel

Simple handling and maintenance



Three levels of filtration – solid, foam and HEPA filter

Capacity of collection vessel for day long treatment

Soft wheels with breaks for easy transport between treatment rooms

Technical specification	DO M	
Nominal voltage, frequency	230, 50/ 60	V, Hz
Maximum vacuum	12	kPa
Max. flow capacity	1100	lit.min-1
Operating mode	S1-100%	
Input power	3,1	A
Noise level	≤48	dB(A)
Weight	32	kg
Dimensions (w x l x h)	350x575x890	mm